

### Product Description





















HandyBelt is used for transfer support. Most commonly, the belt is applied on the patient, and the caregiver then uses the belt as an aid in the transfer. Alternatively, it can be applied on both patient and caregiver so that the patient also can use the belt as support when transferring.

#### **IMPORTANT!**

Transferring a patient always involves a certain level of risk. Read the instructions for use before use. It is important to completely understand the contents of the instructions for use. The equipment should only be used by trained personnel. Ensure that the accessory is suitable for the transfer. Exercise care and caution during use. As a caregiver, you are always responsible for the patient's safety. You must be aware of the patient's ability to make it through the transfer. If something is unclear, contact the manufacturer or supplier.

## Symbol Description

These symbols can be found in this document and/or on the product.

Symbol	Description
	Warning. Used where extra care and attention is needed.
	Read instructions for use before use
	CE mark
	Legal manufacturer
	Date of manufacturing
	Product Identifier
	Serial Number
	Medical device
	Do not dry clean
	Do not iron
	Tumble drying
	Do not tumble dry
	Washing Symbol is always in combination with recommended temperature in °C and °F.
	Latex Free
	Periodic inspection
	Do not bleach
	Line drying in shade
	Direction upwards and symbol on outside of equipment.
	Symbol on outside of equipment.
 <small>(01) 0100887761997127 (11) YYMMDD (21) 012345678910</small>	GS1 Data Matrix Barcode that may contain following information. (01) Global Trade Item Number (11) Production Date (21) Serial Number

## Safety Instructions

### Before use, keep the following points in mind:

- Plan the transfer operation so that it can be done as safely and smoothly as possible.
- A responsible person within your organization should decide on a case-by-case basis whether one or more caregivers are needed.
- The patient should understand and follow the instructions so he/she can be guided to a natural pattern of movement during raising.
- The patient should have head and torso stability and be able to sit upright. When gait training, it is also important that the patient has good balance when standing and the ability to move his/her feet independently when walking.
- Before transferring, check that the HandyBelt is applied correctly.
- HandyBelt is not a lifting aid and the handles are not intended for lifting purposes.
- To minimize the strain for the caregiver, in case the patient should fall, it is important that the caregiver does not by his/her own efforts try to keep the patient in a standing position. Instead, we recommend that the caregiver, in a smooth motion, lower the patient to the floor. Remember to let go of the handles to avoid injuries to the caregiver's hands caused by crushing.
- Always work ergonomically. Use the patient's ability to actively participate in the operations.

 Medical Device Class I Product

### Product Changes

Liko's products undergo continuous development, which is why we reserve the right to make product changes without prior notice. Contact your Liko representative for advice and information about product upgrades.

### Design and Quality by Liko in Sweden

Liko is quality certified according to ISO 9001 and its equivalence for the medical device industry, ISO 13485. Liko is also certified according to environmental standard ISO 14001.

### Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Care and Maintenance

Check the product before each use. Check the following points with regard to wear and damage:

- Fabric
- Straps
- Seams
- Handles
- Safety buckle

### Do not use damaged products.

If anything is unclear, please contact the manufacturer or supplier.

**Washing instructions:** Please see the product label.

### Expected Life Time

Per these Instructions for Use, this product should be inspected regularly, especially after it has been washed. Paying particular attention to wear and damage to seams, straps/handles, fabric, and/or fabric properties. Use in accordance with these instructions and following the instructions for Care and Maintenance, this product will have a continuous useful life so long as the product passes inspection.


### Recycling Instructions

The entire product should be sorted as Combustible waste.


Hill-Rom evaluates and provides guidance to its users on the safe handling and disposal of its devices to aid in the prevention of injury, including, but not limited to: cuts, punctures of the skin, abrasions, and any required cleaning and disinfection of the medical device after use and prior to its disposal. Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories.

If in doubt, the user of the device shall first contact Hill-Rom Technical Support for guidance on safe disposal protocols.

## Application






## Removal




## Areas of use

Raising from sitting to standing, transfer from sitting to sitting and support for standing/walking.



## Overview HandyBelt



Name	Prod. No.	Size	Size Marking (A)	Number of Handles (B)	Length (C/D)
HandyBelt	3760014	S Small	Orange	7	60/100 cm (24"/39")
HandyBelt	3760015	M Medium	Yellow	9	75/120 cm (30"/47")
HandyBelt	3760016	L Large	Blue	11	90/140 cm (35"/55")
HandyBelt	3760017	XL Extra Large	Black	13	110/180 cm (43"/71")

**Material:** polyester

Fabric in Handybelt is made of non-flame retardant materials.

**Maximum load:** Each handle can withstand load up to 30 kg / 66 lbs.

## Alternative Products

It is important to ensure that the transfer operation is not too strenuous for the caregiver. Should the patient need more support, Liko recommends other products from our extensive range of raising and walking aids.

Alternative products in the Liko range:

Sabina sit-to-stand lift

Liko LiftPants in combination with a patient lift

Liko MasterVest in combination with a patient lift



[www.hillrom.com](http://www.hillrom.com)

Liko AB  
Nedre vägen 100  
975 92 Luleå, Sweden  
+46 (0)920 474700

Liko AB is a subsidiary of Hill-Rom Holdings inc.

Enhancing outcomes for  
patients and their caregivers:

**Hill-Rom**