# Hillrom NP50 Monodensity High Resilience Foam Mattress

## Instructions for Use

### Safety and Usage Tips

#### Intended Use

The NP50 mattress benefits are pressure ulcers prevention help. It is suitable for low to moderate risk adult patients, within a recommended weight range of 30 kg and 150 kg.

#### Environment

It can be used as a mattress in the following environments, as defined in the standard IEC 60601-2-52:

- Application environment 1 (acute care);
- Application environment 2 (short-term care in hospitals or other medical establishments);
- Application environment 3 (long-term care in medical establishments);
- Application environment 5 (outpatients or ambulatory care).

#### Contraindications

The NP50 mattress is not designed to come into direct contact with damaged skin and it is not suited to patients with atypical anatomies.

#### Intended Users

The NP series mattresses are designed to be used by Qualified Staff for patient care from several care application environments.

#### First Use

Hill-Rom recommends that the mattress be disinfected before first use.

### User Instructions

**Carefully read all the instructions. Failure to do so could result in equipment damage and injury.**

#### Understanding the Symbols

	Manufacturer
<b>REF</b>	Device Reference
<b>LOT</b>	Batch Code
	CC (Country Code): Country of manufacture Date of Manufacture (YYYY-MM-DD)

<b>MD</b>	Medical Device
<b>UDI</b>	Unique Device Identification
 = 250 kg	Safe Working Load is the greatest permissible load that can be applied after which a damage may occur on the mattress. Safe working load, including the total weight of the patient, accessories and the load supported by these accessories.
<b>CE</b>	Medical Device conformity mark
	Biological Hazard Warning
	<b>WARNING</b> This symbol indicates that the failure to follow the associated recommendation can put the patient or the user in danger, or damage the device.
	Refer to the Instructions For Use
	Does not contain natural latex in accordance with ISO 7000

#### Before Placing the Patient on the Mattress

Remove the NP50 mattress from its polyethylene packaging.

Take care not to damage the mattress with scissors or cutters. Check that the zip fastener is fully closed.

Store the mattress at ambient temperature (18 to 20°C) for at least 24 hours.

Make sure that the mattress is correctly installed and well centered on the sleep surface.

Cover the mattress with a sheet to avoid direct contact between the patient's skin and the surface of the mattress.

#### Complying with Conditions for Use

Symbol	Feature	Use
	Temperature	+10°C to +40°C
	Hygrometry	30% - 85%
	Atmospheric pressure	700 mbar - 1060 mbar

Supervision is advisable if the top edge of the siderail is less than 22 cm above the mattress. When the patient is left unattended, it is recommended to place the bed in the lowest position.

Failure to properly assess the fit of the NP50 mattress and bed frames could result in the risk of entrapment which may cause serious injury or death.

If necessary, place a plate between the patient and the mattress when performing cardio-pulmonary reanimation (external heart massage).

## Transferring the Mattress from one Bed to Another

**Never leave the patient on the mattress during transfer.**

Make sure that the bed chassis brakes are applied to prevent any accidental movements.

If the mattress has handles, use them to transfer the mattress more easily.

While handling the mattress, all necessary precautions must be taken to avoid any damage.

Do not work on the device when it is occupied.

## Cleaning and Disinfecting

Cleaning and disinfection are two separate procedures.

### Complying with Safety Instructions

- Check that the brakes are applied on the bed on which the mattress is installed.
- Lock out all electric functions.
- Never clean the mattress by pouring water on it, nor with a high-pressure jet.
- Never use water at a temperature of more than 70°C.
- Refer to the recommendations of the cleaning and disinfecting of the product provided by the manufacturer.
- Thoroughly dry the mattress before reusing it.
- Suitable individual protective equipment must be worn during the steps of the cleaning operation (blouse, gloves, eye protection, etc.).

**Failure to implement one or more of these recommendations may lead to damage or deterioration, preventing use of the device and rendering the warranty void.**

### Controlling Infection

**Insufficient Cleaning = Risk of Infection (biological danger)!**

All parts must be kept clean at all times in order to avoid the risk of infection. All necessary precautions must be taken to eliminate all visible soiling.

**The following recommendations are not designed to replace cleaning and disinfection protocols drawn up by the hygienic officer or by other bodies for your hospital in the event of particular infectious situations.**

### Complying with Manufacturer Recommendations

The cleaning and disinfecting method described below applies specifically to the mattress and its accessories and is designed to save time and to help combat nosocomial infection more effectively.

**Hill-Rom recommends that the mattress be disinfected before first use.**

**When cleaning, always check the cover for cuts, tears, cracks or snags; never use a mattress with a damaged cover.**

**If the internal foam of the mattress is soiled and presents a risk of infection, then replace the mattress.**

**Following other cleaning and disinfection protocols or using other products than those recommended by Hill-Rom may compromise devices and patient safety, and void the warranty.**

## De-commissioning

Hill-Rom recommends that the mattress should be disinfected before de-commissioning.

Customers should adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories. If in doubt, the user of the device shall first contact Hill-Rom Technical Support for guidance on safe disposal protocols.

### Products To Be Avoided

Never use cleaning agents, detergents, degreasing agents or industrial solvents containing any of the following products in order to avoid damaging the mattress:

	Phenol		Cresol
	Hydrochloric, Nitric or Sulfuric Acids		Soda
	Dimethylformamide		Tetrahydrofuran

Do not use highly acidic detergents or disinfectants (pH<4).

Do not use highly basic detergents or disinfectants (pH>10).

Never use abrasive cleaning materials or products, such as scouring pads.

## Cleaning and Disinfecting after the Departure of a Patient or when Changing the Sheets

### Using Recommended Products

#### Terralin Protect (Schülke)<sup>1</sup>

Concentration: 5 ml/l (= 0.5%), contact time: 1 h.

#### Cleaning

Clean the mattress with a cloth slightly dampened with hot water and a neutral detergent solution. Check that the solution does not contain any of the products to be avoided mentioned above.

#### Cleaning Tough Stains

Quickly wipe away any traces of pharmaceutical solutions used for the patients to avoid permanent damage to the surface.

- Remove tough stains with a soft brush and neutral detergents or a chlorine-based solution according to labeling.
- To eliminate hardened stains (excreta, other forms of soiling), soften them by soaking and take care to thoroughly dry the cover before putting it back on the mattress.

1. **Terralin Protect** is listed by the V.A.H (Verbund für Angewandte Hygiene e.V.).

## Disinfecting

In the event of visible soiling, Hill-Rom recommends that the mattress should be disinfected with an intermediate disinfectant (tuberculocidal) that complies with the applicable regulations (e.g. with the requirements of regulation 2017/745/EEC).

For all other disinfectants:

 $C \leq 1000 \text{ ppm}$	Chlorine-based solutions must be less than or equal to 1,000 ppm. Chlorine-based solutions until 10,000 ppm can be used but will increase discolouration effect.
 25% C <sub>2</sub> H <sub>5</sub> OH	Ethanol (alcohol) can be used. Concentration must not exceed ¼ of ethanol for ¾ of water.

Proceed with the final steps before reusing the cleaned and disinfected mattress (See "Performing the Final Steps", page 3).

## Cleaning and Disinfecting at Regular Intervals or in the Event of High Risk of Contamination

Follow the same instructions as above Go to "Cleaning and Disinfecting after the Departure of a Patient or when Changing the Sheets" on page -2., but with the products listed below.

### Using Recommended Products

#### After the departure of patients with infectious diseases:

- **Dismozon Pur (Bode)**<sup>1</sup>  
 Concentration: 40 g/l (= 4%), contact time: 1 h

#### After the departure of patients with Clostridium Difficile:


- **Dismozon Pur**  
 Concentration: 15 g/l (= 1.5%), contact time: 2 h
- **Anioxy-Spray WS (Anios)**  
 Ready to use, contact time: 30 min

### Dry Steam Cleaning the mattress

Dry steam or superheated steam contains no more than 6% of water in suspension and avoids the effects of condensation.

To avoid damage due to high pressure or an abnormal surface temperature, take the following precautions:

- Do not use accessories such as high pressure hoses (A). It is preferable to use soft non-metallic brushes (B) and wipers (C) in such a way as to reduce the pressure to an acceptable level.






- Do not brush and apply reduced pressure to labels and markings.
- Dry carefully and look for signs of water ingress before reusing.

#### Hill-Rom recommends the use of the Sanivap method:

- Application time recommended by the supplier: approximately 1 sec.
- Distance (recommended by the supplier) between the source and the surface to be disinfected: approximately 3 mm.
- Application: exterior parts of the mattress (covers).
- Frequency: for each new patient.

## Machine Washing the Mattress Cover

Do not always machine-wash, as doing so will reduce the life time of the components. Machine-wash only if particular risks of infection occur.

### Reinstalling the cover

When reassembling the cover after machine washing, make sure to match with the associated foam Go to "Mattress / foam / cover combinations" on page -3.

Orientate properly the cover on the foam regarding the marking located on the foam in the feet area toward the top.

#### Mattress / foam / cover combinations

Ref. Mattress	Ref. Foam	Ref. Cover	Dim. / cm
ASS027	UPH027	TEX027	14 x 85 x 198
ASS028	UPH028	TEX028	14 x 90 x 198
ASS049	UPH049	TEX049	14 x 90 x 203
ASS007	UPH007	TEX007	14 x 85 x 198

### Performing the Final Steps

- Always remove all traces of the products used when washing or disinfecting the mattress.
- There may be a risk of discolouration when using chlorine-based products. Devices must be thoroughly rinsed with a cloth soaked in water, then dried.
- Make sure that all the parts of the mattress are perfectly dry before installing to avoid any risk of condensation forming inside the mattress.

### Complying with Conditions for Transport and Storage

Symbol	Feature	Transport / Storage*
	Temperature	-30°C to +50°C
	Hygrometry	20% - 85%
	Atmospheric pressure	700 mbar - 1060 mbar

\*Applicable only if the mattress is transported or stored in its original packaging.

The mattress must be stored in its original packaging:

- Protected against light and damp;
- At least 10 cm above floor level to prevent fluid ingress;
- Protected against dust;
- Outside passageways.

Never stack more than 15 mattresses.

1. **Dismozon Pur** is listed by the V.A.H (Verbund für Angewandte Hygiene e.V.).

## Specifications

Ref. Mattress	Country	Dim. / cm
ASS027	Europe except UK, IT	14 x 85 x 198
ASS028	Europe except UK, IT	14 x 90 x 198
ASS049	Europe except UK, IT	14 x 90 x 203
ASS007	UK - IT	14 x 85 x 198

The NP50 mattress must only be used with the covers mentioned in the Spare Parts list (see page 4).

Mattress weight: 10 kg.

## Spare Parts

For any information about spare parts, please give the reference printed on the side of the mattress when you call.

## Cover references

Ref. Cover	Type of seam	Ref. Mattress	Dim. / cm
TEX027	Sewn	ASS027	14 x 85 x 198
TEX028	Sewn	ASS028	14 x 90 x 198
TEX049	Sewn	ASS049	14 x 90 x 203
TEX007	Sewn	ASS007	14 x 85 x 198

### NOTE :

*The color of the cover may vary and is non-contractual.*

## Accessories

Part Number	Description
ASS043	NP Mattress Extension (14 x 9 x 85 cm)
ASS044	NP Mattress Extension (14 x 18 x 85 cm)
ASS077	NP Mattress Extension
ASS078	NP Mattress Extension (20 x 21 x 85 cm)

The Mattress Extensions are provided to extend the mattress length at the foot end.

## Regulatory Requirements

 The CE conformity mark applicable to class I Medical Devices was applied for the first time in 2008.

**The mattresses meet the following standards:**

### All mattresses

- EN ISO 10993-1: 2020
- EN ISO 10993-5: 2009
- EN ISO 10993-10: 2010
- EN 597-1:2015
- EN 597-2:2015
- GPEM/CP - Guidelines D1-90, February 1991 (France)

### In addition for the Mattress ASS007

- UNI 9175: 2010
- BS 6807:2006 Clause 9 (Crib5)
- BS 7177:2008 for medium hazard

### Country of manufacture (ISO 3166-1)

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## Contacts

Notice to Users and/or Patients:

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

For any request or to order copies of this manual contact your national Hill-Rom representative or go to [hillrom.com](http://hillrom.com) and order the article with the part number 163791.

Hillrom is a trademark of Baxter International Inc. or its subsidiaries.  
Terralin is a registered trademark of Schülke & Mayr GmbH.  
Dismozon is a registered trademark of Bode Chemie GmbH.  
Sanivap is a registered trademark of Sanivap S.A..

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Hillrom er et varemerke for Baxter International Inc. eller dets datterselskaper.
Terralin er et registrert varemerke for Schülke & Mayr GmbH.
Dismozon er et registrert varemerke for Bode Chemie GmbH.
Sanivap er et registrert varemerke for Sanivap S.A.
Hill-Rom S.A.S.
Hill-Rom S.A.S. oppfyller ISO 13485- og ISO 9001-standardene