**Welch Allyn® Elite Retinoscope 18240, 18245, 18300**Instructions for use  
771381 DIR 8002597 Rev. B Revision date: 2020-08

REF # 901024 RETINOSCOPE

**Rx ONLY**

© 2020 Welch Allyn, Inc. All rights reserved. No one is permitted to reproduce or duplicate in any form, this instructions for use or any part thereof without permission from Welch Allyn. Welch Allyn® is a trademark of Welch Allyn, Inc. Hillrom™ is a trademark of Hill-Rom Services, Inc.

Welch Allyn, Inc. ("Welch Allyn") assumes no responsibility for any injury to anyone that may result from (i) failure to properly use the product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual, or (ii) any illegal or improper use of the product.

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
hillrom.com

Welch Allyn, Inc. is a subsidiary of Hill-Rom Holdings, Inc.

EC REP and EU IMPORTER  
Welch Allyn Limited  
Navan Business Park, Dublin Road  
Navan, Co. Meath C15 AW22  
IrelandAuthorized Australian Sponsor  
Welch Allyn Australia Pty Ltd.  
Unit 401, 2-4 Lyndon Park Road  
Macquarie Park, NSW 2113  
Phone 1800 650 083**ENGLISH****Instructions for use**

Thank you for purchasing the Welch Allyn® 3.5V Retinoscope. This instrument has been designed to meet the needs of today's practitioners and incorporates features not found on any other retinoscopy:

- External Focusing Sleeve - Unique planetary gear system allows for easy adjustment no matter what size hand or how instrument is held. Continuous 360° rotations. Maintains the same plane of focus during rotation.
- Improved Light Output - Brighter halogen lamp provides 50 percent more intensity than previous lamps. The reflex is now crisper and easier to see in all patient's eyes. Retinoscopy can be done faster and more accurately.
- Dust-free Optics - New housings and glass cover on the front keep the instrument cleaner longer.
- Crossed Linear Polarizing Filter - Dramatically reduces glare from lenses. Allows retinoscopy to be performed closer to the eye and with confidence.
- Fixation Cards - New cards that easily attach increase the ease with which dynamic retinoscopy is performed.
- Improved Optics - Glare and shadows have been eliminated for a clearer and more precise view.
- Interchangeability - By simply changing the lamp, the streak retinoscopy can be converted to a spot retinoscopy.

**Indications for use**

The Retinoscope is intended to aid in measuring the refraction of the eye. The Retinoscope is intended for use by clinicians and medically qualified personnel. Welch Allyn Retinoscopes are intended for use in both ambulatory and acute care settings.

**Symbol descriptions**For information on the origin of these symbols, see the Welch Allyn symbols glossary: [welchallyn.com/symbols/glossary](http://welchallyn.com/symbols/glossary)

MD Medical device # Product Identifier

Manufacturer EC REP Authorized Representative in the European Community

Lot code Consult Instructions for use

REF Reorder number GTIN Global Trade Item Number

Temperature limitation % Humidity limitation

Rx ONLY Prescription only or For Use by or on the order of a licensed medical professional

Separate collection of Electrical and Electronic Equipment. Do not dispose as unsorted municipal waste.

CAUTION The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

WARNING The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death. Warning symbols will appear with a grey background in a black and white document.

Residual risk This product complies with relevant electro-magnetic interference, mechanical safety, performance, and biocompatibility standards. However, the product cannot completely eliminate potential patient or user harm from the following:

- Harm or device damage associated with electro-magnetic hazards.
- Harm from mechanical hazards.
- Harm from device, function, or parameter unavailability.
- Harm from misuse such as inadequate cleaning, and/or
- Harm from device exposure to biological triggers that may result in a severe systemic allergic reaction.

**Product use**

The operation of the control sleeve of the scope The operator will note that the width of the streak varies as the sleeve is raised and lowered (see Figure 1). When the operating sleeve is in the lowest position the light rays emitted are slightly divergent. The instrument acts with a plano mirror effect, which reflects divergent rays that will never come to a focus. As the sleeve is raised, the streak focuses. With the sleeve all the way up, the retinoscope acts with a concave mirror effect, where the light rays cross and then diverge. Below the rays cross, the eye's reflex moves in opposite directions with the concave mirror effect as compared to the plano mirror effect.

Throughout this IFU, we will use the plano mirror effect unless specified.

The rotary movement of the control sleeve mechanism allows the streak to rotate 360° to ascertain the axis of astigmatism (Figure 1).

**Preliminary Steps**

- Set the sleeve in its lowest position (plano-mirror effect).
- Position yourself 213 centimeters (2 meters) from the patient. This distance implies a working lens of +1.50D (computed as the reciprocal of working distance in meters). Working distance and lens may be varied to suit the practitioner's needs. (In this instruction book the 2/3 meter (26') working distance is assumed. Different working distances can be used, but remember to adjust for your working distance.)
- With the refracting equipment in place, direct the practitioner's attention to a fixation spot at 15 feet or more from the eye and align the streak vertically as the streak or in the opposite direction.

- Observe the "reflex" which will appear as in Figure 2, providing no oblique astigmatism present. If oblique astigmatism is present the reflex will appear more like Figure 3, where the reflex does not appear vertical.
- Move the vertical streak horizontally across the pupil and back again and observe whether the reflex moves in the same direction as the streak or in the opposite direction.


Figure 2

Figure 3

6. Rotate the control sleeve until the streak is horizontal and move the streak vertically. The reflex will appear as in Figures 4 or 5.

Figure 4

Figure 5

7. If the streak and the reflex move in the same direction with no lens in the refractive apparatus, the refraction is one of these:

- Hyperopia;

- Emmetropia;
  - Myopia of less than 1.50 diopters.
- If the reflex moves in the opposite direction, the error is myopia greater than 1.50 diopters.

**Determining refractive error by neutralization**

Before starting, make sure the eye not being refracted has some "against" motion using the plano mirror effect. This will blur vision to prevent accommodation. If "with" or neutral motion is noticed initially, place about a +0.00 sphere before the eye once neutral motion is seen.

**Neutralizing with spheres only**

- Change shape in the minus direction until the reflexes in all axes have "with" motion.
- Adjust in the plus direction until the reflex fills the pupil in one meridian and all motion has stopped. This will be one of the principal meridians if astigmatism is present. That meridian is then said to be neutralized.

**Test for neutralization by one of these methods:**

- Move the sleeve all the way up (concave mirror position); the reflex should also appear neutralized;
- Move closer to the patient and "with" motion should return; or
- Place an extra +0.25 sphere in the apparatus and "against" motion should appear;

4. Repeat the neutralization in the meridian 90° away.

**Locating the axis of astigmatism**

Two phenomena in determining the axis of astigmatism: break and drift. Break is observed when the streak is not aligned with a principal meridian of the astigmatism (Figures 3 and 5). The streak will be aligned with a principal meridian when the break effect disappears and the width of the reflex is narrowest (and it appears its brightest, Figure 6).



Figure 6  
Proceed with neutralization as before-neutralizing one principal meridian first, then 90° away to neutralize the second principal meridian (Figures 7 and 8).

Figure 7


Figure 8

**Interpretation of Results**

- Hyperopia
  - Hyperopia exists when, at the 213-meter distance using the plano mirror effect, "with" motion is neutralized using a plus lens greater than +1.50 diopters and both meridians neutralize with the same strength lens.
  - Total hyperopia is estimated by subtracting 1.50 diopters from the total strength lens used. For example, if it takes a +2.50 lens to neutralize motion at 213 meter, the total hyperopic error is +1.00 diopter.

- Myopia (exists under several circumstances)
  - When neutralizing, using the plano mirror effect at 2/3 meter, the neutralization will be less than 1.50 diopter strength. (When motion is neutralized with exactly a 1.50 diopter lens, the eye is emmetropic.)
  - When at 2/3 meter, using the plano mirror effect, no motion appears at all. In other words, when the motion is neutralized with no lens in the refracting apparatus. The myopia is then exactly 1.50 diopters.
  - When the motion is "against" using the plano mirror effect, and is neutralized with a minus lens.

- Astigmatism
  - Astigmatism exists when the two principal meridians neutralize with different strength lenses. It may be present in many forms.
  - Simple hyperopic;
  - Compound hyperopic;
  - Compound myopic;
  - In the mixed form (one meridian hyperopic and the opposite one myopic).

- Astigmatism can be measured in one of two ways:
  - Neutralize one principal meridian first. Then add the appropriate plus or minus cylindrical lens until the other principal meridian is also neutralized by continuing to add spherical lenses until the second principal meridian is neutralized. Then the astigmatic error is equal to the difference in strength of lenses necessary to neutralize the two meridians.
  - Neutralize the cylinder by continuing to add spherical lenses until the second principal meridian is neutralized. Then the astigmatic error is equal to the difference in strength of lenses necessary to neutralize the two meridians.

**Special Considerations**

Axis of astigmatism: Extreme care must be used in setting the axis of the cylinder, if the correcting cylinder is of the proper power, a 10° error in axis will produce a new astigmatism of approximately one third of the strength of the original astigmatism with its principal meridian at approximately 45° to those of the original astigmatism. The technique for setting the axis is referred to as "straddling". When you have an approximate correction of the refractive error and wish to refine the axis setting, the following technique will be helpful:

- Move up closer to the eye so that the edge of the reflex can be seen, and compare the widths of the two reflexes as you rotate the streak 45° to either side of the correcting cylinder axis. Recede slowly while doing this.

- Compare the widths of the two reflexes. If there is an axis error, the reflex will be different in the two positions. If you are using plus cylinders, rotate the axis toward the narrow band until the reflex widths are equal.

- With minus cylinders, move the axis away from the narrow band. When the reflex widths are equal, the proper axis has been determined. It is important that the spherical and cylindrical strength be checked again after completion of this maneuver.

**Spiralling**

Spiralling is a method for estimating astigmatism without lenses. This can be helpful in determining the starting point for lens introduction when working with patients who have a high unknown astigmatism (Figure 9).

The Welch Allyn Streak retinoscope is particularly well suited for this technique because the instrument maintains the same focal plane during streak rotation.

**Crossed linear polarizing filter (CLPF)**

This filter can be engaged by moving the sliding switch on the practitioner's side of the instrument from the down to the up position (Figure 9).

This filter cuts down on reflections and allows retinoscopy to be performed closer to the axis of the correcting lens.

**Instructions for lamp replacement**

No. 18240 and 18245 Streak Retinoscopes (Use only Hill Rom 3.5V Halogen Lamp No. 08200-U)

18300 Spot Retinoscope (Use only Hill Rom 3.5V Halogen Lamp No. 08300-U)

1. Remove retinoscope from power source (Figure 10).


Figure 10

2. Remove lamp: Lift out with nail file, letter opener, or similar instrument under base flange (Figure 11).


Figure 11

3. Insert new lamp:
 

- 08200-U lamp: line up pin on lamp with slot between the electrical contact wires (Figure 12); lamp base contact receptacle as far as it will go.
- 08300-U lamp: push lamp into receptacle as far as it will go. Lamp should insert easily; do not force. Lamp base contact pin should be even with metal cutouts in retinoscope base (Figure 13).

4. Replace retinoscope on power source.

• The Streak and Spot retinoscopes are essentially the same instrument. The No. 18240 &amp; 18245 Streak retinoscopes can be converted to a Spot retinoscope by simply inserting the 08300 lamp, and vice versa.

**Cleaning, disinfection, and maintenance**

• RetinoScope external housings may be cleaned with a mild detergent and soft cloth.

• Windows may be cleaned with alcohol on a cotton swab or lens paper.

• Hyperopia;

• Emmetropia;

• Myopia of less than 1.50 diopters.

If the reflex moves in the same direction with no lens in the refractive apparatus, the refraction is one of these:

- Hyperopia;

**WARNING Do not immerse.**

- Fixation cards may be cleaned with a mild detergent.
- Routinely inspect the device for wear and damage. Do not use if you see signs of damage, if the instrument malfunctions, appears not to be working properly, or if you notice a change in performance. Contact Hillrom Technical Support for guidance on safe disposal protocols.

**Disposal**

Users must adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories. If in doubt, the user of the device shall first contact Hillrom Technical Support for guidance on safe disposal protocols.

**PATENT/PATENTS**

hillrom.com/patents

May be covered by one or more patents. See above Internet address. The Hill-Rom companies are the proprietors of European, US, and other patents and pending patent applications.

**Hillrom Technical Support**For information about any Hillrom product, contact Hillrom Technical Support at [hillrom.com/en-us/about-us/locations](http://hillrom.com/en-us/about-us/locations).**Notice to users and/or patients in the EU**

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Standards**

This product complies with ISO 12865, ISO 15004-1, ISO 10993-1, IEC 60601-1 and IEC 60601-1-2.

Country-specific standards are included in the applicable Declaration of Conformity.

**Lot code**

Lot code YY-JJJ where YY = last 2 digits of the year, JJJ= consecutive day of the Julian year.

**Accessories**

18250 Fixation card set

08200-U 3.5V Halogen HPX TM Streak lamp

08300-U 3.5V Halogen HPX TM Streak retinoscope lamp

**Warranty**

Welch Allyn warrants that these retinoscopes will be free of defects in workmanship or for one year from the date of purchase (excludes lamps and batteries).

Welch Allyn will correct any of these defects at its expense by repairing or replacing the product. This warranty does not cover damage caused by accident, repairs performed by a provider who is not authorized by Welch Allyn, or failure to follow Welch Allyn's maintenance instructions.

THIS LIMITED WARRANTY IS PROVIDED INSTEAD OF ANY OTHER IMPLIED PRODUCT QUALITY WARRANTIES AND IT IS THE SOLE WARRANTY OF THE PRODUCTS

Figure 2

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