## Welch Allyn Neonatal 1-5 Blood Pressure Cuff

Instructions for use

**[REF]** 80026303 Ver. B. Revised 2020-04

These instructions apply to  [#] 901045 blood pressure cuff, single patient.

<span></span>	Consult Instructions for use <p>Consulter le mode d'emploi</p> <p>Consulte las instrucciones de uso</p> <p>Gebrauchsanweisung beachten</p> <p>Consultare le istruzioni per l'uso</p> <p>Beoordeling de gebruiksaanwijzing</p> <p>Se bruksanvisningen</p> <p>Der hvenvis til brugsanvisningen</p> <p>Tutustu käyttöohjeeseen</p> <p>Läs bruksanvisningen</p> <p>Ανατρέξτε στις οδηγίες χρήσης</p> <p>Consultar as instruções de utilização</p> <p>Καταγράψτε τα στοιχεία βασικών</p> <p>См. руководство по эксплуатации</p> <p>Zapoznać się z instrukcją obsługi</p>	<span></span>	One Port <p>Un port</p> <p>Ein puerto</p> <p>Ein Anschluss</p> <p>Uma porta</p> <p> Een poort</p> <p>En port</p> <p>En port</p> <p>Yksi portti</p> <p>En port</p> <p>Μία θύρα</p> <p>Uma conexão</p> <p>Tek Port</p> <p>Один порт</p> <p>Jeden port</p>	<span></span>	Two Ports <p>Deux ports</p> <p>Zwei puertos</p> <p>Zwei Anschlüsse</p> <p>Due porte</p> <p>Twee poorten</p> <p>To porter</p> <p>To porte</p> <p>Kaksi porttia</p> <p>Tuä portar</p> <p>Δύο θύρες</p> <p>Dois conexões</p> <p>Tek Port</p> <p>Два порта</p> <p>Dua porty</p>
---------------	---	---------------	---	---------------	---

<span></span>	Use larger cuff <p>Utiliser un brassard plus grand</p> <p>Utilice un manguito más grande</p> <p>Größere Manschette verwenden</p> <p>Utilizzare un bracciale più grande</p> <p>Groteret manchet gebruiken</p> <p>Bruk en større mansett</p> <p>Brug større manchet</p> <p>Käytä suurempaa mansettia</p> <p>Använd större manschett</p> <p>Χρησιμποιήστε μεγαλύτερο χρομό</p> <p>Utilizati braccialele maiori</p> <p>Geniş kaf kullanın</p> <p>Используйте манжету побольше</p> <p>Użyj większy mankiet</p>	<span></span>	Global Trade Item Number <p>Référence de commerce international</p> <p>Número de artículo de comercio global</p> <p>GTIN (Global Trade Item Number)</p> <p>Número artículo per il commercio globale</p> <p>Global Trade Item Number</p> <p>Globalt handelsartikelnummer</p> <p>Globalt varennummer</p> <p>Kansainvälinen tuotenumero</p> <p>GTIN-artikkelnummer</p> <p>Γlobaλitítraktiartiklónumero</p> <p>Número de Artigo Comercial Global</p> <p>Global Ticari Ürün Numarası</p> <p>Глобалны нумар товару</p> <p>Globally unique handover</p>	<span></span>	Limb Circumference (Min./Max) <p>Circonférence du membre (min./max.)</p> <p>Circunferența de la extremitate (Min./Măx.)</p> <p>Umfang der Extremität (min./max.)</p> <p>Circonférence de l'artère (min./max.)</p> <p>Omtrek van arm/ben (min./max.)</p> <p>Lemtomérés (min./maks)</p> <p>Kroppsdelens omkreds (min./maks)</p> <p>Raajan ympärämitta (minimi/maksimi)</p> <p>Kroppsdelens omkrets (min./max)</p> <p>Περιφέρεια του μέμβρου (Min./Máx.)</p> <p>Use Circumf. (min./maks.)</p> <p>Окружность конечности (минималная/максималная)</p> <p>Obwód kończyny (min./maks.)</p>	<span></span>	Not made with natural latex <p>Non fabriqué à partir de latex naturel</p> <p>Non fabbricato con latex naturale</p> <p>Nicht aus Naturlatex hergestellt</p> <p>Non realizzato in lattice naturale</p> <p>Ikke fremstillet af naturligt latex</p> <p>Innehåller ikke naturligt latex</p> <p>Indeholder ikke naturligt latex</p> <p>Etä sisältä luonnollista lateksiä</p> <p>Innehåller inte naturligt latex</p> <p>Δεν έχει κατασκευαστεί από φυσικό λάτεξ</p> <p>Não fabricado com látex natural</p> <p>Dogal lateleke yapılmamıştır</p> <p>Не содержит натурального латекса</p> <p>Wykonano bez użycia naturalnego lateksu</p>
---------------	---	---------------	--	---------------	---	---------------	--

<span></span>	Product identifier <p>Identifiant du produit</p> <p>Identificador de producto</p> <p>Prod-ID</p> <p>Identificativo del prodotto</p> <p>Productificatie</p> <p>Prod-ID</p> <p>Produkt-ID</p> <p>Tuotetunnus</p> <p>Productidentifikation</p> <p>Αυτοματισμένο προϊόν</p> <p>Identificador do produto</p> <p>Ürün kimliği</p> <p>Идентификатор изделия</p> <p>Identyfikator produktu</p>	<span></span>	Index Edge <p>Bord du repère</p> <p>Borde del índice</p> <p>Indekstante</p> <p>Bordo índice</p> <p>Indeksgrens</p> <p>Indeksleyerpunkt</p> <p>Indekstant</p> <p>Merkin reun</p> <p>Indekstant</p> <p>Espoz óbrazit</p> <p>Extremidade do índice</p> <p>Indeks Eşjiği</p> <p>Край указателя</p> <p>Krawędź wskaźnika</p>	<span></span>	Artery Index <p>Repère d'indice artériel</p> <p>Markador de índice arterial</p> <p>Arterienindexmarkierung</p> <p>Indicatore dell'arteria</p> <p>Slagaderindexmarkering</p> <p>Arterieleksmärker</p> <p>Arterie-Indexmarke</p> <p>Vältnorm kohdistusmerkki</p> <p>Arterindeksmärki</p> <p>Δείκτης αρτηριακού οπίσθιου</p> <p>Markador indicador de arteria</p> <p>Arter indeks (jarekyleyiti)</p> <p>Мерка-указателя артерии</p> <p>Znacznik wskaźnika tętnicy</p>
---------------	--	---------------	---	---------------	--

<span></span>	Reorder number <p>Número de commande</p> <p>N° de reposition</p> <p>Nachbestellnummer</p> <p>Numero per il riordino</p> <p>Bestelnummer</p> <p>Beställingsnummer</p> <p>Beställingsnummer</p> <p>Tilbæsningsnummer</p> <p>Beställingsnummer</p> <p>Αριθμός επαναπληρωτικής παραγγελίας</p> <p>Número de encomenda</p> <p>Yeniiden sipariş numarası</p> <p>Номер повторного заказа</p> <p>Numar повторноgo zamowienia</p> <p>nia</p>	<span></span>	Medical Device <p>Dispositif médical</p> <p>Producto sanitario</p> <p>Medizinprodukt</p> <p>Dispositivo medic</p> <p>Medisch hulpmiddel</p> <p>Medisinsk erhet</p> <p>Medisinsk udstyr</p> <p>Läkemedletslinjen laite</p> <p>Medicinteknisk produkt</p> <p>Λογισμικό υγειονομικού ουσιαστί</p> <p>Dispositivo médico</p> <p>Tiðbiti Chaz</p> <p>Медицинское устройство</p> <p>Wyrob medyczny</p>	<span></span>	Range <p>Intervallo</p> <p>Bereich</p> <p>Intervallo</p> <p>Bereik</p> <p>Område</p> <p>Område</p> <p>Alttavall</p> <p>Ευρος</p> <p>Intervallo</p> <p>Intervallo</p> <p>Diapazon</p> <p>Zakres</p>
---------------	---	---------------	--	---------------	--

<span></span>	Stacking limit by number <p>Limite d'empilement par chiffre</p> <p>Limite de aplado por número</p> <p>Maximale Stapelhöhe nach Anzahl</p> <p>Limiti di accostamento per numero</p> <p>Stapelgrens op aantal</p> <p>Stabilitätsbegrensning etter nummer</p> <p>Stabilitätsbegrensning efter nummer</p> <p>Määränsäilytysrajat</p> <p>Staplingsgräns, antal</p> <p>Όριο σταθμής από αριθμό</p> <p>Limite de empilhamento por número</p> <p>Sayıyla üst üste koyma sın</p> <p>Допустимая высота штабелирования (в шт.)</p> <p>Limit ustawienia opakowań jedno na drugim</p>	<span></span>	Humidity limitation <p>Limite d'humidité</p> <p>Limite de humedad</p> <p>Zulässige Luftfeuchtigkeit</p> <p>Limite di umidità</p> <p>Lochtigheidsbeëik</p> <p>Luffuktigheitsgrens</p> <p>Grenser for luffuktighet</p> <p>Sallittu ilmastotilasto</p> <p>Luffuktigheitsgrens</p> <p>Όριο υγρασίας</p> <p>Limites de humidade</p> <p>Nem szemlélt</p> <p>Диапазон влажності</p> <p>Dopuszczalna wilgotność</p>
---------------	--	---------------	---

<span></span>	Temperature Limit <p>Limites de température</p> <p>Limite de temperatura</p> <p>Zulässiger Temperaturbereich</p> <p>Limite di temperatura</p> <p>Temperaturlgrens</p> <p>Temperaturlgrens</p> <p>Temperaturlgrens</p> <p>Temperaturlgrens</p> <p>Temperaturlgrens</p> <p>Όριο θερμοκρασίας</p> <p>Limites de temperatura</p> <p>Скалікi Smm</p> <p>Температурний діапазон</p> <p>Dopuszczalna temperatura</p>	<span></span>	Manufacturer <p>Fabricante</p> <p>Hersteller</p> <p>Produttore</p> <p>Fabrikant</p> <p>Produzent</p> <p>Producent</p> <p>Valmistaja</p> <p>Tilvikare</p> <p>Κατασκευαστής</p> <p>Fabricante</p> <p>Uretici</p> <p>Производитель</p> <p>Produsător</p>	<span></span>	LOT <p>Lot code number</p> <p>Número de code du lot</p> <p>Numero de código de lote</p> <p>Chargencode-Nummer</p> <p>Lot code</p> <p>Número de codice lotto</p> <p>Lotkodennummer</p> <p>Lotkodennummer</p> <p>Eden tunnusnumero</p> <p>Αριθμός κωδικού παρτίδας</p> <p>Número de código de lote</p> <p>Lot Kodu Numarası</p> <p>Номер кода партии</p> <p>Numar kodu partii</p>
---------------	---	---------------	---	---------------	---

<span></span>	Authorized Representative in the European Community <p>Représentant autorisé dans la Communauté européenne</p> <p>Representante autorizată en la Comunitatd Europea</p> <p>Autorisierter Händler in der Europäischen Gemeinschaft</p> <p>Rappresentante autorizzato nella Comunità europea</p> <p>Geautoriseerde vertegenwoordiger in de Europese Gemeenschap</p> <p>Autorisert representant i EU</p> <p>Goðkendt representant í EÚ</p> <p>Les autorizovaný predstaviteľ v Európskom spoločenstve</p> <p>Autorizowany przedstawiciel we Wspólnocie Europejskiej</p>	<span></span>	Single patient - multiple use <p>Patient unique - utilisation multi-tiple</p> <p>Usò múltiple en un solo paciente</p> <p>mehrfache Verwendung für nur einen Patienten</p> <p>Singolo paziente - Uso multiplo</p> <p>Meeroudig gebruik bij één patiënt</p> <p>Kun til bruk på én pasient</p> <p>Enkelt patient - flergangsbrug</p> <p>Yhden potilaan toistuva käyttö</p> <p>En patient - flera användningar</p> <p>Ενας ασθενής - πολλαπλή χρήση</p> <p>Um patiente - varias utilizações</p> <p>Tek hastada kullanan i'cinder</p> <p>Для использования только с одним пациентом</p> <p>Jeden pacjent – wielokrotny użytek</p>
---------------	---	---------------	--

For information on the origin of these symbols, see the Welch Allyn symbols glossary: welchallyn.com/symbolglossary

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician or licensed health-care practitioner. This device should be used by trained personnel.

## ENGLISH

### Introduction

**Intended use**
Welch Allyn Neonatal Blood Pressure Cuffs are noninvasive blood pressure cuffs intended for use in conjunction with automated sphygmomanometers capable of neonatal measurement to determine blood pressure in neonatal patients. Welch Allyn Blood Pressure Cuff is intended to be used by a clinically trained, medical professional.

**Contraindications**
Welch Allyn Neonatal Blood Pressure cuffs are contraindicated for use with non-automated sphygmomanometers or automated sphygmomanometers without a neonatal mode.

**Residual risks**
This product complies with relevant electro-magnetic interference, mechanical safety, performance, and biocompatibility standards. However, the product cannot completely eliminate potential patient or user harm from the following:

- Harm or device damage associated with electro-magnetic hazards
- Harm from mechanical hazards
- Harm from device, function, or parameter unavailability
- Harm from misuse error, such as inadequate cleaning, and/or
- Harm from device exposure to biological triggers that may result in a severe systemic allergic reaction

**Compatible devices**
This blood pressure cuff is designed to function with automated sphygmomanometers with neonatal mode tested to:
ANSI/AAMI/ISO 81060-2:2013, Non-invasive sphygmomanometers — Part 2: Clinical investigation of automated measurement type.

### Warnings

A warning statement in this manual identifies a condition or practice which, if not corrected or discontinued immediately could lead to injury, illness, or death.

**WARNING:** Do not reuse a cuff on more than one patient. Reusing a cuff on more than one patient might cause crosscontamination.

**WARNING:** Do not apply cuff to areas on patient where skin is delicate or damaged. Check cuff site frequently for irritation.

**WARNING:** Possible measurement error. Only use the correct size cuff applied as directed and using included/appropriate connectors in conjunction with automated equipment capable of neonatal measurement.

**WARNING:** Do not over tighten cuff.

**WARNING:** Do not apply cuff to limbs used for IV infusion.

**WARNING:** Patient injury risk. Never install luer lock connectors on Welch Allyn blood pressure tubing.

Using these connectors on blood pressure cuff tubing creates the risk of mistakenly connecting this tubing to a patient's intravenous line and introducing air into the patient's circulatory system.

**WARNING:** Minimize cuff movement and limb motion during readings.

**WARNING:** Intravenous Systems (IV) - Do not connect cuff's connectors to any intravenous fluid systems.

**WARNING:** Only use cuff with automated sphygmomanometers capable of neonatal mode BP measurement.

### Cautions

A caution statement in this manual identifies conditions or practices that could result in damage to the equipment or other property.

### Operation

Use the Welch Allyn Neonate Blood Pressure cuff the same as a traditional blood pressure cuff. The blood pressure cuff works with automated sphygmomanometers that include a neonatal measurement mode.

Operational Pressure Range: 0 to 150 mmHg (Neonatal mode)

- Select the cuff size appropriate for the patient's arm circumference. The applicable range in centimeters is printed on each cuff. The artery index marker must fall within the Range indicated on the cuff.
- Apply the cuff to the patient's limb, verify the size, and close, ensuring a snug but not a tight or loose fit.
  - If the artery index marker falls short of the range, use a larger cuff to ensure accurate results.
  - If the artery index marker falls beyond the range, use a smaller cuff to ensure accurate results.
- Connect the cuff to an automated sphygmomanometer.

### Cleaning and low-level disinfection

**WARNING:** Cleaning and disinfection procedures must be conducted by persons trained in medical device cleaning and disinfection.

**WARNING:** Consult the cleaning and germicidal cleaner agents' manufacturer's instructions for their proper use and germicidal efficacy.

**CAUTION:** Prevent liquid from entering cuff tubing and bladder or damage to the cuff bladder may occur. **CAUTION:** Use only the cleaning or germicidal cleaner agent types listed or damage to cuff may occur. **CAUTION:** Improper reprocessing may cause degradation of device; follow inspection procedures to assure integrity of device.

**CAUTION:** Do not aggressively scrub cuff; damage to cuff markings and/or cuff closure integrity may occur.

### Materials

- Neutral pH enzymatic-based cleaning detergent.
- Bleach-based germicidal cleaner suitable for use on healthcare equipment and capable of low-level disinfection. For example: A cleaner including 1:10 dilution of bleach (6500 ppm sodium hypochlorite) and detergent. Reference EPA-registered disinfectants: http://www.epa.gov/oppad001/chemregindex.htm.
- Clean or sterile cloths, soft brush, soaking tray, and potable rinse water (softened preferred).

### Cleaning only (cuffs and connector accessories)

- Prepare neutral pH enzymatic cleaning detergent solution per manufacturer instructions.
- Immerse or soak cuff and accessories in solution. Do not allow liquid ingress into bladder tubing.
- Soft brush all surfaces of the cuff and accessories to remove visible soil. Repeat as necessary.

### Clean and disinfect (cuffs)

- Clean: Thoroughly saturate (spray or immerse) all surfaces of the cuff and accessories with germicidal cleaner.
  - Soft brush to remove visible soil.
  - Water rinse. Do not allow liquid ingress into bladder tubing.
  - Camp dry and inspect.
- Disinfect: Thoroughly re-saturate (spray or immerse) all surfaces of the cuff and accessories with germicidal cleaner.
- Soft brush all surfaces. Allow a 5 minute wet contact time or longer, if directed by the germicidal cleaner manufacturer. Do not exceed 10 minutes of wet contact time.

### After cleaning or cleaning and disinfection

- Thoroughly water rinse. Do not allow liquid ingress into bladder tubing.
- Damp dry with a clean cloth
- Allow to air dry.
- Inspect cuff for deterioration, adequate closure integrity, and inflate to assess for leaks. Do not use if any abnormalities are found.

### Technical support

For information about any Welch Allyn product, contact Hillrom Technical Support: hillrom.com/en-us/about-us/locations

**Lot code (marked directly on cuff)**
**Lot decoder:** YYDDT where
YY= Last two digits of the year
DDD= Day of the year
T= Machine

## Environmental specifications

NOTE: Comply with regional law when non-automated sphygmomanometer or accessories are disced	
Storage temperature	-20°C to 70°C (-4°F to 158°F)
Operating temperature	10°C to 40°C (50°F to 104°F)
Operating relative humidity	15% to 90% (noncondensing)

### Safe disposal

Users must adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices and accessories. If, in doubt, the user of the device shall first contact Hillrom Technical Support for guidance on safe disposal protocols

### Patent/patents

May be covered by one or more patents. See below Internet address. The Hill-Rom companies are the proprietors of European, US, and other patents and pending patent applications.

hillrom.com/patents

### Legal statement

Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal and improper use of the product, that may result from failure to properly use the product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

### Notice to users and/or patients in EU

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Warranty

Your Welch Allyn™ product, when new, is warranted to be free from original defects in material and workmanship and to perform in accordance with manufacturer's specifications under normal use and service.

© 2020 Welch Allyn, Inc. All rights reserved. No one is permitted to reproduce or duplicate, in any form, this instructions for use or any part thereof without permission from Welch Allyn. Welch Allyn™ and FlexiPort™ are trademarks of Welch Allyn, Inc. Hillrom™ is a trademark of Hill-Rom Services, Inc.

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