## Hillrom.

### Welch Allyn® FlexiPort® EcoCuff® Blood pressure cuff for single patient use

<b>Part number</b>	901045 Blood Pressure Cuff, Single Patient
<b>REF</b>	771352, 80026304 Ver. B
<b>Revision date:</b>	2020-04

**Caution** US federal law restricts this device to sale by or on the order of a physician or licensed healthcare practitioner. This device should be used by patients only.

**Standards**
This blood pressure cuff is designed to function within the limits prescribed by:
ANSI/AAMI/ISO 81060-1:2007 (ENISO 81060-1:2012)
ENISO 81060-1:2012
Non-invasive sphygmomanometers—Part 1: Requirements and test methods for non-automated measurement type.
2.
Non-invasive sphygmomanometers—Part 1: Requirements and test methods for non-automated measurement type.

#### Hillrom technical support

hillrom.com/en-us/about-us/locations

© 2020 Welch Allyn. All rights are reserved. To support the intended use of the product described in this publication, the purchaser of the product is permitted to copy this publication, for internal distribution only, from the media provided by Welch Allyn. No other use, reproduction, or distribution of this publication, or any part of it, is permitted without written permission from Welch Allyn. Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

Welch Allyn®, FlexiPort®, and EcoCuff® are trademarks of Welch Allyn, Inc. Hillrom™ is a trademark of Hill-Rom Company, Inc. For further information, visit hillrom.com/patients.

### ENGLISH

#### Introduction

The FlexiPort EcoCuff blood pressure cuff is intended to be used on the upper arm in conjunction with noninvasive blood pressure measurement systems. The cuff is nonsterile and is available in pediatric through adult sizes. The device is not intended for neonatal use. The cuff and the FlexiPort EcoCuff blood pressure cuffs is not designed, nor intended for use, except as indicated. Welch Allyn Blood Pressure Cuff is intended to be used by a clinically trained medical professional.

**Contraindications**
Welch Allyn EcoCuff blood pressure cuffs are contraindicated for neonate use. Do not use on neonate patients.

**Compatible devices**
This blood pressure cuff is designed to function with nonautomated sphygmomanometers tested to: ANSI/AAMI/ISO 81060-1:2007 (ENISO 81060-1:2012)
ENISO 81060-1:2012
Non-invasive sphygmomanometers - Part 1
Requirements and test methods for nonautomated measurement type.
This blood pressure cuff is designed to function with automated sphygmomanometers tested to: ANSI/AAMI/ISO 81060-1:2013.
Non-invasive sphygmomanometers - Part 2
Clinical investigation of automated measurement type and automatication.

#### Warnings

The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death.

**WARNING** Intravenous (IV) systems—Do not connect cuffs with luer lock connectors to IV fluid systems. Otherwise, air might enter the patient's blood stream.

**WARNING** Do not reuse a cuff on more than one patient. Reusing a cuff on more than one patient might cause cross-contamination.

**WARNING** Do not clean, disinfect, or sterilize the cuff. These methods have not been tested for effectiveness or impact to cuff performance.

**WARNING** Use the cuff only when the port falls within the range indicator window on the cuff. Otherwise, erroneous readings might result.

**WARNING** Do not apply the cuff to areas on the patient where skin is delicate or damaged. Check the cuff site frequently for irritation.

**WARNING** Do not apply the cuff to an arm used for IV infusion. Inflating the cuff on this arm might interrupt IV flow.

**WARNING** Minimize cuff movement and arm motion during readings. Excessive movement might affect readings.

**WARNING** Ensure an airtight seal at all connection points before use. Excessive leaking might affect readings.

**WARNING** When applying the cuff to the patient, allow space for 1 to 2 fingers between the patient's upper arm and the cuff to ensure a proper fit.

**WARNING** Keep debris out of the port on the cuff. Debris in the cuff can damage the cuff.

**WARNING** Intravenous (IV) systems—Do not connect cuffs with luer lock connectors to intravenous fluid systems. Otherwise, fluid might enter the cuff and cause infection.

**Caution** The caution statements in this manual identify conditions or practices that could damage the equipment or other property.

**Caution** Close the hook-and-loop fastener before inflating the cuff. Failure to properly close the fastener might lead to cuff damage.

**Residual risk**
This product complies with the relevant electro-magnetic interference, mechanical safety, performance, and biocompatibility standards. However, the product cannot completely eliminate potential or user harm from the following:

- Harm from mechanical hazards.
- Harm from device use when a parameter is unavailable.
- Harm from misuse error, such as inadequate cleaning, and/or
- Harm from device exposure to biological triggers that may result in a severe systemic allergic reaction.

**Operation**
Notes:
• The cuff works with one-tube and two-tube manual and automated sphygmomanometers.
• The cuff works only with FlexiPort fittings.
• Use the cuff on one or more times on one patient and then discard it.

- Do not reuse the cuff on another patient.
- Do not attempt to clean, sterilize, or disinfect the cuff.

To take a blood pressure measurement:

- To apply the cuff, align the artery marker **1** on the cuff with the brachial artery and wrap the cuff around the upper arm, so that the part of the cuff with the range indicator window **2** is on top of the port **3**.

- Ensure that the port **4** shows through the range indicator window **2** and connect the FlexiPort fitting **5** to the port through the range indicator window. (If the port does not fall within the range indicator window, select a different cuff size.)

- To remove the cuff, disconnect the FlexiPort fitting from the port.
- Unwrap the cuff from the upper arm.

**Do not mark directly on cuff!**

YY-MM-DD TT
Where Y = Year, M = Month, D = Day, and T = Tool

**Environmental specifications**
**Caution** Store and use the cuff within stated temperature and humidity ranges. Otherwise, the cuff might not meet performance specifications.

Storage temperature: -20°C to 70°C (4°F to 158°F)
Storage relative humidity: 15% to 95% (noncondensing)
Operating temperature: 10°C to 40°C (50°F to 104°F)
Operating relative humidity: 15% to 90% (noncondensing)

**Safe disposal**
Users must adhere to all federal, state, regional, and/or local laws and regulations as it pertains to the safe disposal of medical devices.
**WARNING** Dispose of this device in a manner that meets the user of the device's local technical support for guidance on safe disposal protocols.

**Notice to Users and/or Patients in EU:**
This device and its accessories are not to be used unless the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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