Welch Allyn® PanOptic™ Ophthalmoscope

Welch Allyn® Ophtalmoscope PanOptic™

Welch Allyn® Oftalmoscopio PanOptic™

118 Series / Série 118 / Serie 118
(VET)

Instructions for use

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REF 773582, DIR 80028104 Ver. A, Revision date: 2021-03

©2021 Welch Allyn, Inc. All rights are reserved.
Printed in the U.S.A.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
hillrom.com

Welch Allyn, Inc. is a subsidiary of Hill-Rom Holdings, Inc.

Authorized Australian Sponsor
Welch Allyn Australia Pty. Ltd.
Unit 4.01, 2-4 Lyonpark Road
Macquarie Park, NSW 2113
Phone 1800 650 083

ENGLISH

Symbols

For information on the origin of these symbols, see the Welch Allyn symbols glossary: welchallyn.com/symbolsglossary

WARNING The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death. Warning symbols will appear with a grey background in a black and white document.

CAUTION The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

NOTE Presents clarification about an instruction or helpful information about a feature or behavior.

Manufacturer

TYPE BF – Indicates this is a product with Type BF applied parts (the patient eyecup).

Warnings and cautions

WARNING The light emitted from the ophthalmoscope is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument, when operated at maximum intensity, will exceed the safety guideline after 13 minutes when lit with a LED lamp (blue base). Exposure times are cumulative for a 24-hour period.

WARNING When using Halogen illumination, no acute optical radiation hazards are identified. However, Welch Allyn recommends limiting the intensity of the light directed into the patient's eye to the minimum level that is necessary for diagnosis. Infants, aphakes, and persons with diseased eyes are at a greater risk. The risk is also increased if the examined person has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source within the previous 24 hours. This will apply particularly if the eye has had exposure to retinal photography. The intended use of this device is for routine ophthalmic exams on the order of typically less than 60 seconds per eye. Although there is a benefit versus risk factor in any medical procedure, these more complicated exams should not exceed a three minute exam time in 24 hours. Significant use of this device beyond its intended use is not recommended; it may cause harm to the eyes.

WARNING To minimize lamp housing temperatures, on-time should not exceed 2 minutes with off-time not less than 10 minutes.

WARNING Use only PanOptic Patient Eyecups with this product.

WARNING Rx only (for professional use).

WARNING Lamps may be hot. Lamp should be allowed to cool before removal.

WARNING Equipment not suitable for use in the presence of a flammable anesthetic.

CAUTION Halogen lamps are pressurized to provide maximum efficiency and illumination. Mistreatment may cause shattering. Protect lamp surface against abrasion and scratches. Verify power is off when replacing lamp. Dispose of lamp with care.


NOTE A graph of the relative spectral output of the instrument is available upon request.

The PanOptic™ Ophthalmoscope

Transparency of the cornea, lens, and vitreous humor permits the practitioner to directly view arteries, veins, and the optic nerve of the retina.

Direct observation of the structures of the fundus through a PanOptic Ophthalmoscope may show disease of the eye itself or may reveal abnormalities indicative of disease elsewhere in the body. Among the most important of these are vascular changes due to diabetes or hypertension and swelling of the optic nerve head due to papilledema or optic neuritis. In this sense, the eye serves as a window through which many valuable clinical evaluations may be made.

When a preliminary diagnosis of an imminently dangerous eye condition, such as acute (angle-closure) glaucoma or retinal detachment, is made by the examiner, prompt referral to an eye specialist may prevent irreversible damage. Or, when distressing but less urgent conditions, such as visual impairment due to cataract or vitreous floaters, are recognized, the patient can be reassured and referred.


A	Patient Eyecup	D	Aperture/Filter Dial
B	Patient's Side	E	Focusing Wheel
C	Soft Grip Handle	F	Practitioner's Side Brow Rest

Set-up instructions

Attach the PanOptic Ophthalmoscope head to a Welch Allyn power source

- Align cutouts in the PanOptic Ophthalmoscope base with lugs on power source. (The PanOptic Ophthalmoscope fits all 3.5v Welch Allyn power sources.)
- Push in slightly and turn 90 degrees.


Attach the Patient Eyecup to the PanOptic

WARNING Use only PanOptic Patient Eyecups (part number 118092) with this product.

- Insert eyecup into the patient's side of the instrument.
- Push and twist in one motion until you feel the two "lock" together.

See Accessory and Replacement Parts List to reorder new PanOptic Patient Eyecups.

Attach the Corneal Viewing Lens to the PanOptic (model 11820 only)

- Insert Corneal Viewing Lens into patient's side of the instrument.
- Push and twist in one motion until you feel the two "lock" together as shown, except replacing Patient Eyecup with the Corneal Viewing Lens.
- Verify that the line on the lens housing is aligned with the top seam of the PanOptic instrument.

Apertures and filters

There is a wide range of practical apertures and filters to select from: small spot, large spot, micro spot, slit aperture, red-free filter, cobalt blue filter (optional), and half-moon aperture (optional).


- Small aperture:** Provides easy view of the fundus through an undilated pupil. Always start the examination with this aperture and proceed to micro aperture if pupil is particularly small and/or sensitive to light. This position is the "Home" position on the aperture dial and is denoted by the green marking.
- Large aperture:** Standard aperture for dilated pupil examination of the eye.
- Micro spot aperture:** Allows easy entry in very small, undilated pupils.
- Slit aperture (models 11810, 11820 only):** Helpful in determining various elevations of lesions, particularly tumors and edematous discs.
- Red-free filter (model 11810, 11820 only):** This filter excludes red rays from the examination field: this is superior to ordinary light in viewing slight alterations in vessels, minute retinal hemorrhages, ill-defined exudates and obscure changes in the macula. The nerve fibers become visible and the observer may note the disappearance of such fibers, as in optic nerve atrophy. The background appears gray, the disc appears white, the macula appears yellow, the fundus reflex is intense and the vessels appear almost black. This filter is also used to help distinguish veins from arteries; veins stay relatively blue, but oxygenated arterial blood makes arteries appear blacker. This greater contrast makes differentiation easier for the examiner.
- Cobalt blue filter (model 11820 only):** In conjunction with fluorescein dye applied topically to the cornea and an add-on magnifying lens (included), this filter is helpful in detecting corneal abrasions and foreign bodies. In this way, it can be used as a surrogate to a Woods Light.
- Half-moon aperture (11810 only):** Provides a combination of depth perception and field-of view.

Lamp replacement instructions

WARNING Lamps may be hot. Lamp should be allowed to cool before removal. Use only Welch Allyn model 03800 lamp with this product.

CAUTION Halogen lamps are pressurized to provide maximum efficiency and illumination. Mistreatment may cause shattering. Protect lamp surface against abrasion and scratches. Verify power is off when replacing lamp. Dispose of lamp with care..

To insert the new lamp, align tab on lamp with notch in the base of the ophthalmoscope head and push inward until firmly seated.

Cleaning and repair instructions

Cleaning and repairs (with exception of the lamp replacement) beyond what is specified in the "general cleaning instructions" section should be performed at the factory.

General cleaning instructions

NOTE Excess solution entering the optical assembly could damage internal components. Use caution to ensure cloth is not saturated with solution.

This is a precision optical instrument. Please follow these specific cleaning requirements:

- Do not use any solvent based cleaners on the lenses. Use only a clean/soft optical grade cleaning cloth on the lenses.
- The PanOptic Ophthalmoscope is Non-ster ilizable.
- Disinfect housing using soft cloth moistened with alcohol.
- Clean the Patient Eyecup with soft cloth/cotton swab moistened with alcohol.
 - Soft cloth/cotton swab moistened with alcohol.
 - Soaking in gluteraldehyde based solutions (i.e. Cidex OPA).

NOTE For disinfection time / procedure requirements it is recommended that you follow the directions of the specific brand you are using.

- Autoclaving.
- Do not use acetone based products or other harsh chemicals.


No other preventative maintenance is necessary for this product.

Troubleshooting

Trouble area	Possible cause	Corrective action
No light output	Aperture dial is in-between positions	Rotate the aperture dial.
	Lamp burned out.	Replace the lamp using Welch Allyn part number 03800.
	Wrong lamp/ incorrect Welch Allyn lamp/non-Welch Allyn lamp.	Replace lamp using Welch Allyn part number 03800.
Spot is not centered.	The battery handle has completely discharged.	Charge the handle, check the battery, and/or the charger.
	The aperture dial is not centered.	Move aperture dial to the full detent position.
Not achieving full or expected view.	Unit is not in proper operating position.	Ensure that the eyecup is slightly compressed during the procedure.
Cannot obtain sharp focus/hazy view.	The lenses are dirty.	Clean the lenses with a clean/soft optical grade cleaning cloth.

Trouble area	Possible cause	Corrective action
Dim light output	Lamp pin is not engaged in slot	Insert the lamp with pin firmly seated in slot.
	Film on optics	Clean the lenses with a clean/soft optical grade cleaning cloth.
	The handle is not fully charged.	Charge the handle, check the battery and/or the charger.
Glare in view	Fingerprints or dirt on the objective (patient's side) lens	Clean the objective (patient's side) lens with a clean/soft optical grade cleaning cloth.
Cannot attach to power handle	Lamp is not fully inserted.	Insert the lamp with pin firmly seated in slot.

Accessory and replacement parts list

Item	Description	Product No.
1	Corneal Viewing Lens (model no. 11820 only)	11875
2	Patient Eyecups (Package of five (5) each 118092)	11870
3	Patient's Side Bumper	118051
4	Lamp	03800
5	Chrome Ring	118027
6	Practitioner's Side Brow Rest	118052

Specifications

Model number

11810 PanOptic Ophthalmoscope without Cobalt Blue Filter
11820 PanOptic Ophthalmoscope with Cobalt Blue Filter and Add-on Corneal Viewing Lens

Dimensions

5.125"L x 1.4"W x 3.750"H without eyecup
Eyecup 1.45"L

Weight

.48 lbs without eyecup
.50 lb with eyecup

Equipment classification

IPX0 Equipment not protected against the ingress of water.

WARNING Equipment not suitable for use in the presence of a flammable anesthetic.

Environment

Transport/Storage:-

-20 °C - 49 °C, 95% Rh max, 500hPa - 1060 hPa, altitude

Operating:

10 °C - 49 °C, 95% Rh max, 500 hPa - 1060 hPa, altitude

Standards and compliance

EN60601-1
IEC 60601-1-2
CAN/CSA-C22.2 No. 601.1-M90
UL 2601-1, Second Edition, 1997

Warranty

PanOptic Ophthalmoscope Warranty

Welch Allyn, Inc. warrants the No. 11810, 11820 PanOptic Ophthalmoscope to be free of original defects in material and workmanship and to perform in accordance with manufacturer's specifications for a period of one year from the date of purchase. If this instrument or any component thereof is found to be defective or at variance with the manufacturer's specifications during the warranty period, Welch Allyn will repair or replace the instrument or component(s) at no cost to the purchaser.

This warranty only applies to instruments purchased new from Welch Allyn or its authorized distributors or representatives. The purchaser must return the instrument directly to Welch Allyn or an authorized distributor or representative and bear the costs of shipping.

This warranty does not cover breakage or failure due to tampering, misuse, neglect, accidents, modification or shipping, and is void if the instrument is not used in accordance with manufacturer's recommendations or if repaired or serviced by other than Welch Allyn or a Welch Allyn authorized representative.

No other express or implied warranty is given.

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