



Hillrom™

IMPROVING MOVEMENT AND INDEPENDENCE



Your guide to finding the right
Patient Handling Solutions for
your clients' needs.



[CLICK HERE TO VIEW FULL SCREEN](#)





Hillrom™

WELCOME TO THE HILLROM SAFE PATIENT HANDLING EQUIPMENT E-BOOK

We understand you want to specify the best and most appropriate equipment for your clients, but the healthcare landscape and patients' needs are complex and the solutions are many.

To help you navigate which lift, sling, vest or accessory might be the most appropriate to complete the task you and your client require, we have created this E-Book.

It's quick and easy to use, simply click the relevant answer to a set of questions and we will guide you to the most appropriate equipment required to complete your task.

For the final stage, simply add your contact details and one of our Hillrom representatives will be in touch to discuss in further detail the perfect solution for your clients' needs.

[CLICK HERE TO GET STARTED](#)



CONTENTS



05	Plus size solutions Hillrom Ultra RepoSheet™ Hillrom UltraSling™ Hillrom Ultra LiftPants™
-----------	---

20	Turning / repositioning with a reduced care package Hillrom RepoSheet™ Hillrom HandySheet™
-----------	---

21	Turning / repositioning with multiple carers Hillrom RepoSheet™ Hillrom HandySheet™
-----------	--


22-26	Lifting Solutions with a reduced care package Hillrom HygieneVest™ Hillrom HygieneSling™ Hillrom Original HighBack Sling™ Hillrom SilhouetteSling™ Hillrom ComfortSling™ Plus High Hillrom Soft Original HighBack Sling™ Hillrom UniversalSling™
--------------	--

27-31	Lifting Solutions with multiple carers Hillrom HygieneVest™ Hillrom HygieneSling™ Hillrom Original HighBack Sling™ Hillrom SilhouetteSling™ Hillrom ComfortSling™ Plus High Hillrom Soft Original HighBack Sling™ Hillrom UniversalSling™
--------------	---

32-33	Standing transfers Hillrom SupportVest™ Hillrom SafetyVest™ Hillrom ComfortVest™
--------------	--

34	Standing & stepping transfers with a reduced care package Hillrom LiftPants™ Hillrom MasterVest™
-----------	---

35	Standing & stepping transfers with multiple carers Hillrom LiftPants™ Hillrom MasterVest™
-----------	--





Hillrom™

Do you need a solution for a plus size client?



Yes



No





Hillrom™

Plus size solutions



Hillrom Ultra RepoSheet™ with a Hillrom UltraTwin Overhead Lift

OVERVIEW

- Facilitates repositioning a client higher up in bed and turning in bed with a lift, with minimal burden on the caregiver.
- Can be used in combination with or instead of sheets.
- Intended for use with clients who cannot actively participate in repositioning in bed, or alternatively, are too heavy for the caregivers to reposition.

PRODUCT DETAILS

Regular: 255 x 135cm

Ultra: 255 x 145cm

Fabrics: 50% Cotton / 50% Polyester, Net Polyester

Fabrics: 50% Cotton / 50% Polyester, Net Polyester

SWL: 200kg

SWL: 500kg



Can be used in combination with the Hillrom UltraTwin™ Overhead lift system

Hillrom UltraSling™ with a Hillrom UltraTwin Overhead Lift

OVERVIEW

- Developed for lifting plus sized people in a safe and comfortable way.
- Supports the entire back up to the shoulders.
- Leg supports are connected individually on the sling bar, which provides extra space for the abdomen and legs, for optimal support.

PRODUCT DETAILS

SWL: 500kg

Fabric: Polyester



Can be used in combination with the Hillrom UltraTwin™ Overhead lift system

Hillrom Ultra LiftPants™ with a Hillrom UltraTwin Overhead Lift

OVERVIEW

- Supports the whole body and gives heavier clients the confidence to take those first steps.
- For clients who can put weight on their legs. Allows natural pattern of movement when rising.
- Safe and secure for raising, standing, gait and balance training.
- Takes the strain in the crotch area and by the hip loop straps.

PRODUCT DETAILS

SWL: 500kg

Fabrics: Net Polyester

Sizes Available: XXL, XXXL



Can be used in combination with the Hillrom UltraTwin™ Overhead lift system

[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)





Hillrom™

Reposition / turn in bed | Full lift sitting transfer | Stand & / or stepping



What task are you trying to complete with your client?



Reposition / turn in bed



Full lift sitting transfer



Stand & / or stepping





Hillrom™

Is there a reduced care package or multiple carers completing the task?



Reduced care package



Multiple carers



Hillrom™

Is there a reduced care package or multiple carers completing the task?



Reduced care package



Multiple carers





Hillrom™


Is there a reduced care package or multiple carers completing the task?



Reduced care package



Multiple carers





Hillrom™

Standing transfer | Standing & stepping transfer



**Will your client
be completing
a standing
transfer?
Or a standing and
stepping transfer?**



Standing transfer



Standing & stepping transfer






Hillrom™

Full lifting seating transfer with reduced package



Does your client need to undress during this task?



Undressing the client



Not undressing the client






Hillrom™

Full lifting seating transfer with multiple carers 


Does your client need to undress during this task?



Undressing the client



Not undressing the client





Hillrom™

Full seating transfer with a reduced care package - undressing



Does your client need head support from the sling for this task?



Head support needed for client



No head support needed for client





Hillrom™

Full seating transfer with a multiple carers - undressing



Does your client need head support from the sling for this task?



Head support needed for client



No head support needed for client





Hillrom™

Full lifting with a reduced care package in a wet area / sling left in situ



Does your client need head support from the sling for this task?



Head support needed for client



No head support needed for client





Hillrom™

Full lifting with a multiple carers in a wet area / sling left in situ



Does your client need head support from the sling for this task?



Head support needed for client



No head support needed for client





Will the client's arms be inside or outside of the vest?



Outside the vest



Inside the vest





Hillrom™

Full seating transfer with a reduced care package - not undressing 


Will the sling be removed after the transfer, or is a breathable fabric required to be left in-situ / for wet environments?



Sling will be removed after transfer is complete



Sling should be breathable for left in-situ or wet environments





Hillrom™

Will the sling be removed after the transfer, or is a breathable fabric required to be left in-situ / for wet environments?



Full seating transfer with multiple carers - not undressing



Sling will be removed after transfer is complete



Sling should be breathable for left in-situ or wet environments





Hillrom™

Turning and / or repositioning with a reduced care package



Hillrom RepoSheet™ with a Hillrom Overhead Lift

OVERVIEW

- For clients who cannot actively participate in repositioning in bed, or are too heavy for the caregivers to reposition.
- Can be used in combination with or instead of bedsheets.
- Works best when used in combination with an overhead lift.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Regular: SWL 200kg

Fabrics: 50% Cotton / 50% Polyester or Net Polyester

Ultra: SWL 500kg

Fabrics: 50% Cotton / 50% Polyester or Net Polyester



Hillrom HandySheet™

OVERVIEW

- **HandySheet™** and **HandyTube™** are manual assistive devices that facilitate transfers by reducing friction between the client and the underlying surface.
- The sliding sheets in the Handy™ series are made of a high-quality, low-friction fabric that slides easily, even after long-term use.

PRODUCT DETAILS

Handytube short:

70cm (W) x 90cm (L) (28 x 35")

Handysheet long, M:

90cm (W) x 185cm (L) (35 x 73")

Handysheet long, XL:

140cm (W) x 215cm (L) (55 x 85")



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Turning and / or repositioning with multiple carers



Hillrom RepoSheet™ with a Hillrom Mobile lift

OVERVIEW

- For clients who cannot actively participate in repositioning in bed, or are too heavy for the caregivers to reposition.
- Can be used in combination with or instead of bedsheets.
- Hillrom mobile lifts have a range of Safe Working Loads and lift heights. The Golvo™ 9000 mobile lift is recommended for use with the Hillrom RepoSheet™ due to its high lift height and flexible liftstrap.

PRODUCT DETAILS

Regular: SWL 200kg

Fabrics: 50% Cotton / 50% Polyester or Net Polyester



Hillrom HandySheet™

OVERVIEW

- **HandySheet™** and **HandyTube™** are manual assistive devices that facilitate transfers by reducing friction between the client and the underlying surface.
- The sliding sheets in the Handy™ series are made of a high-quality, low-friction fabric that slides easily, even after long-term use.

PRODUCT DETAILS

Handytube short:

70cm (W) x 90cm (L) (28 x 35")

Handysheet long, M:

90cm (W) x 185cm (L) (35 x 73")

Handysheet long, XL:

140cm (W) x 215cm (L) (55 x 85")



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with a reduced care package, undressing mid-transfer with head support



Hillrom HygieneVest™ with a Hillrom Overhead Lift

OVERVIEW

- Specifically designed to safely and securely lift people with low muscle tone to the toilet.
- Suitable for people with high / low tone and spasticity.
- Upright sitting posture.
- Gives good support for back and head.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, M, L, XL

SWL: S - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with a reduced care package, undressing mid-transfer with no head support



Hillrom HygieneSling™ with a Hillrom Overhead Lift

OVERVIEW

- A simple and easy-to-use sling which suits most clients. It facilitates dressing / undressing and for transfers to and from the toilet.
- Easy to apply.
- For people with good torso stability.
- The client holds arms outside the sling.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, M, L, XL

SWL: S - L = 200kg, XL = 300kg



Hillrom HygieneSling™ with safety belt with a Hillrom Overhead Lift

OVERVIEW

- **HygieneSling™ Mod. 45** is equipped with a **SafetyBelt™** which lifts around the upper body for extra stability.
- Equipped with **SafetyBelt™**, suitable for clients with poor torso stability.
- Can be used with one of the client's arms on the inside.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, M, L, XL

SWL: S - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)


Hillrom™

 Lifting with a reduced care package in a wet area / sling left in situ with head support
 

Hillrom Original HighBack Sling™ with a Hillrom Overhead Lift

OVERVIEW

- Enables a reclined sitting posture with head support, suitable for clients with poor head and torso stability.
- Arms are positioned inside the sling.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL
Fabric: Net Polyester

SWL: S - L = 200kg, XL = 500kg



Hillrom SilhouetteSling™ with a Hillrom Overhead Lift

OVERVIEW

- Provides a reclined sitting posture with comfortable head support.
- Specifically designed for transfers to narrow wheelchairs where the sling should be left in the chair.
- Thin and pliable without handles or reinforcements.
- Ideal for people with postural deformities.
- Client can rest their arms on the inside or outside.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: XS, S, M, L
Fabric: Net Polyester

SWL: 200kg



Hillrom ComfortSling™ Plus High with a Hillrom Overhead Lift

OVERVIEW

- Supports the entire back with a soft head support.
- Applied and removed with the client in a supine position and remains in the chair after the transfer.
- Ideal for clients with arthritis / leg amputation.
- Possibility of open hip joint angle combined with a 4 point sling bar.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, M, L, XL
Fabric: Net Polyester

SWL: S - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with a reduced care package in a wet area / sling left in situ, with no head support



Hillrom UniversalSling™ with a Hillrom Overhead Lift

OVERVIEW

- Basic model which fits most clients.
- Excellent for older clients with limited mobility as enables an upright sitting posture with support for the entire back.
- Net polyester material makes this sling a feasible option for bathing and other wet environments.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL

SWL: XS - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)





Hillrom™

Lifting with a reduced care package - sling removed after transfer



Hillrom Soft Original HighBack Sling™ with a Hillrom Overhead Lift

OVERVIEW

- Comfortable fit with convenient head support when in a seated position.
- Allows a reclined sitting posture with head support, ideal for clients with poor head and torso stability.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL, XXL

SWL: S - L = 200kg, XL - XXL = 500kg



Hillrom UniversalSling™ with a Hillrom Overhead Lift

OVERVIEW

- Basic model which fits most clients. Enables an upright sitting posture with support for the entire back.
- The client's arms are positioned inside the sling.

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg. The lift motors can be supported by installed rails or by a FreeSpan gantry system.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL

SWL: S - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with multiple carers and undressing mid-transfer with head support



Hillrom HygieneVest™ with a Hillrom M230 Mobile Lift

OVERVIEW

- Intended for safe lifting and transfer of clients to and from the toilet, or when access to the client is required for hygiene purposes and head support is needed.

Hillrom M230 Mobile Lift

- Ideal for nursing homes and home environment. This model is excellent for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.
- The **Hillrom M230 Mobile Lift** has a SWL of 182kg, weighing a total of 40kg and has maximum lift height of 1,725mm.

PRODUCT DETAILS

Sizes: XXS, XS, S, M, L, XL

SWL: XXS - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with multiple carers and undressing mid-transfer with no head support



Hillrom HygieneSling™ with a Hillrom M230 Mobile Lift

OVERVIEW

- A simple and easy-to-use sling which suits most clients. It facilitates dressing / undressing for transfers to and from the toilet.
- Easy to apply.
- For people with good torso stability.
- The client holds arms outside the sling.

PRODUCT DETAILS

Sizes: S, M, L, XL
SWL: S - L = 200kg, XL = 300kg

Hillrom M230 Mobile Lift

- Ideal for nursing homes and home environment. This model is suited for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.
- The **Hillrom M230 Mobile Lift** has a SWL of 182kg, weighing a total of 40kg and has maximum lift height of 1,725mm.



Hillrom HygieneSling™ with safety belt with a Hillrom M230 Mobile Lift

OVERVIEW

- **HygieneSling™ Mod. 45** is equipped with a **SafetyBelt™** which lifts around the upper body for extra stability.
- Equipped with **SafetyBelt™**, suitable for clients with poor torso stability.
- Can be used with one of the client's arms on the inside.

PRODUCT DETAILS

Sizes: S, M, L, XL
SWL: S - L = 200kg, XL = 300kg

Hillrom M230 Mobile Lift

- Ideal for nursing homes and home environment. This model is suited for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.
- The **Hillrom M230 Mobile Lift** has a SWL of 182kg, weighing a total of 40kg and has maximum lift height of 1,725mm.



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with multiple carers in a wet area / sling left in situ with head support



Hillrom Original HighBack Sling™ with a Hillrom M230 Mobile Lift

OVERVIEW

- Provides a reclined sitting posture with head support, an ideal option for clients with poor head and torso stability.
- Suitable for lifting to and from the floor.
- Net polyester material makes this sling a feasible option for bathing.

Hillrom M230 Mobile Lift

Ideal for nursing homes and home environment.

This model is ideal for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL
Fabric: Net Polyester

SWL: S - L = 200kg, XL = 500kg



Hillrom SilhouetteSling™ with a Hillrom M230 Mobile Lift

OVERVIEW

- Provides a reclined sitting posture with comfortable head support.
- Specifically designed for transfers to narrow wheelchairs where the sling should be left in the chair.
- Thin and pliable without handles or reinforcements.
- Net polyester material makes this sling a feasible option for bathing.

Hillrom M230 Mobile Lift

Ideal for nursing homes and home environment.

This model is ideal for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.

PRODUCT DETAILS

Sizes: XS, S, M, L
Fabric: Net Polyester

SWL: 200kg



Hillrom ComfortSling™ Plus High with a Hillrom M230 Mobile Lift

OVERVIEW

- Supports the entire back with a soft head support.
- Applied and removed with the client in a supine position and remains in the chair after the transfer.
- Ideal for clients with arthritis / leg amputation.
- Possibility of open hip joint angle combined with a 4 point sling bar.

Hillrom M230 Mobile Lift

Ideal for nursing homes and home environment.

This model is ideal for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.

PRODUCT DETAILS

Sizes: S, M, L, XL
Fabric: Net Polyester

SWL: S - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with multiple carers in a wet area / sling left in situ with no head support



Hillrom UniversalSling™ with a Hillrom M230 Mobile Lift

OVERVIEW

- Basic model which fits most clients.
- Excellent for older clients with limited mobility as enables an upright sitting posture with support for the entire back.
- Net polyester material makes this sling a feasible option for bathing and other wet environments.

Hillrom M230 Mobile Lift

Ideal for nursing homes and home environment.

This model is ideal for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL

SWL: XS - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Lifting with multiple carers - sling removed after transfer



Hillrom Soft Original HighBack Sling™ with a Hillrom M230 Mobile Lift

OVERVIEW

- Enables a reclined sitting posture with head support, ideal for clients with poor head and torso stability.
- Suitable for lifting to and from the floor.

Hillrom M230 Mobile Lift

Ideal for nursing homes and home environment.

This model is ideal for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL, XXL

SWL: S - L = 200kg, XL -XXL = 500kg



Hillrom UniversalSling™ with a Hillrom M230 Mobile Lift

OVERVIEW

- Excellent for older clients with limited mobility.
- Enables an upright sitting posture with support for the entire back.

Hillrom M230 Mobile Lift

Ideal for nursing homes and home environment.

This model is ideal for daily transfers of adults and children; for instance, for transfers to and from a wheelchair, toilet and floor.

PRODUCT DETAILS

Sizes: S, MS, M, L, XL

SWL: S - L = 200kg, XL = 300kg



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Standing transfer - arms outside sling



Hillrom SupportVest™ with a Hillrom Sabina II EE Sit-to-Stand Lift

OVERVIEW

- A comfortable sit-to-stand vest which is simple to use and requires little individual setting.
- Ideal for clients who need assistance to raise themselves but who can then stand firmly on their legs.
- Lifts behind the back and under the arms.

PRODUCT DETAILS

Sizes: S, M, L, XL
SWL: 200kg



Hillrom Sabina Sit-To-Stand lift

- Designed for clients who have difficulty standing up without assistance from a seated position.
- Ideal for clients who are able to actively participate in the raising motion.
- When standing, clients can be moved to a wheelchair or toilet; this gives standing practice in connection with the transfer.

The **Sabina Sit-to-Stand Lift** has a SWL of 200kg, weighing a total of 45.2kg.




Hillrom SafetyVest™ with a Hillrom Sabina II EE Sit-to-Stand Lift

OVERVIEW

- A sit-to-stand vest which is unique in safety, form and function.
- Lifts behind the back and around the torso.
- The self-tensioning safety belt eliminates the risk of the client sliding out of the vest.

PRODUCT DETAILS

Sizes: M, L
SWL: 200kg



Hillrom Sabina Sit-To-Stand lift

- Designed for clients who have difficulty standing up without assistance from a seated position.
- Ideal for clients who are able to actively participate in the raising motion.
- When standing, clients can be moved to a wheelchair or toilet; this gives standing practice in connection with the transfer.

The **Sabina Sit-to-Stand Lift** has a SWL of 200kg, weighing a total of 45.2kg.



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Standing transfer - arms inside sling



Hillrom ComfortVest™ with a Hillrom Sabina II Comfort Sit-to-Stand Lift

OVERVIEW

- Differs in a decisive way from other sit-to-stand vests.
- The outer vest lifts behind the back and outside the arms, the inner safety belt prevents sliding and therefore contributes to increased safety.

PRODUCT DETAILS

Sizes: S, M, L, XL
SWL: 200kg

Hillrom Sabina Sit-To-Stand lift

- Designed for clients who cannot tolerate any pressure underneath their shoulder.
- Designed for clients who have difficulty standing up without assistance from a seated position.
- Ideal for clients who are able to actively participate in the raising motion.
- When standing, clients can be moved to a wheelchair or toilet; this gives standing practice in connection with the transfer.

The Sabina Sit-To-Stand Lift has a SWL of 200kg, weighing a total of 45.2kg.



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

Standing and stepping transfers with a reduced carer package



Hillrom LiftPants™ with a Hillrom Overhead Lift

OVERVIEW

- With Hillrom LiftPants™ Ambulation Support the client can gait train under secure conditions.
- Provides freedom of movement, while relieving some or all of the burden of body weight.
- Supports safely, allowing the client to move on their own without the risk of falling.

PRODUCT DETAILS

Sizes: S, M, L, XL
SWL: S - L = 200kg, XL = 500kg

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg.

The lift motors can be supported by installed rails or by a FreeSpan gantry system.



Hillrom MasterVest™ with a Hillrom Overhead Lift

OVERVIEW

- Facilitates sit-to-stand movement of clients who are able to support their own weight when standing.
- Ideal aid for standing, gait and balance training.
- When transferring to or from the toilet, the lower body is free for dressing / undressing.
- Option of an additional leg harness to provide additional support through the pelvis.

PRODUCT DETAILS

Sizes: XS, S, M, L, XL, XXL
SWL: XS - L = 200kg, XL - XXL = 300kg

Hillrom Likorall™ Overhead Lift Systems

A range of Safe Working Loads from 200kg to 500kg.

The lift motors can be supported by installed rails or by a FreeSpan gantry system.



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)





Hillrom™

Standing and stepping transfers with multiple carers



Hillrom LiftPants™ with a Hillrom Golvo™

OVERVIEW

- With Hillrom LiftPants™ Ambulation Support the client can gait train under secure conditions.
- Provides freedom of movement, while relieving some or all of the burden of body weight.
- Supports safely, allowing the client to move on their own without the risk of falling.

PRODUCT DETAILS

Sizes: S, M, L, XL
SWL: S - L = 200kg, XL = 500kg

Hillrom Golvo™ 9000

Hillrom mobile lifts have a range of Safe Working Loads and lift heights. The Golvo™ 9000 mobile lift is recommended for use for standing & stepping due to its height, flexible liftstrap and armrests for the client to hold.



Hillrom MasterVest™ with a Hillrom Golvo™

OVERVIEW

- Facilitates sit-to-stand movement of clients who are able to support their own weight when standing.
- Ideal aid for standing, gait and balance training.
- When transferring to or from the toilet, the lower body is free for dressing / undressing.
- Option of an additional leg harness to provide additional support through the pelvis.

PRODUCT DETAILS

Sizes: XS, S, M, L, XL, XXL
SWL: XS - L = 200kg, XL - XXL = 300kg

Hillrom Golvo™ 9000

Hillrom mobile lifts have a range of Safe Working Loads and lift heights. The Golvo™ 9000 mobile lift is recommended for use for standing & stepping due to its height, flexible liftstrap and armrests for the client to hold.



[CLICK HERE TO ARRANGE A CLIENT ASSESSMENT](#)



Hillrom™

ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.

For more information, please contact your local distributor or Hillrom sales representative.

hillrom.com

Hillrom range of lifts and accessories are intended to be used for the lifting and transferring of patients in a variety of care settings. Hillrom slings and accessories are intended to be used in combination with Hillrom lifts for a range of lifting and transferring situations in a variety of care settings.

Class 1

Manufacturer: Liko AB, Nedre vägen 100, 975 92 Luleå, Sweden

These products are regulated health care products which, where required by applicable regulations, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documentation accompanying medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR223601 MAY-2021 ENG – UK

