



Hillrom™



# SAFER. SUPPORTED. SINGLE USE.

**DISPOSABLE PRODUCTS FOR  
HILLROM SAFE PATIENT HANDLING  
& MOBILITY™ SOLUTIONS.**

## Why Choose Single-Patient Use?

Although every patient encounter is unique, we know safety is always your priority. That's why we offer disposable slings, sheets and lift aids as options to help you:

- Reduce the risk of infections spreading<sup>1</sup>
- Meet short-term lifting and transferring needs
- Keep patients comfortable with non-woven and cotton-polyester materials




<sup>1</sup> Infection Control Today. Disposable Products Enhance Infection Prevention Efforts. <https://www.infectioncontroltoday.com/personal-protective-equipment/disposable-products-enhance-infection-prevention-efforts>. Accessed April 29, 2020.



## ORDERING INFORMATION

Part Number	Description	Size	Quantity
3525911-2	Solo HighBack Sling™	M	20 pcs/box
3525915-2	Solo HighBack Sling™	M	20 pcs/box
3525916-2	Solo HighBack Sling™	L	20 pcs/box
3525917-2	Solo HighBack Sling™	XL	8 pc/box
35911915-4	Solo SupportVest™	M	20 pcs/box
35911916-4	Solo SupportVest™	L	20 pcs/box
3683905-4	Solo OctoStretch™ LiftSheet	M	40 pcs/box
3683906 -4	Solo OctoStretch™ LiftSheet	L	40 pcs/box
3695915-2	Solo MultiStrap™	Standard	20 pcs/box
3695916-2	Solo MultiStrap™	Wide	20 pcs/box
36871004	Solo RepoSheet™	Generous	8 pcs/box
36871003-2	Solo RepoSheet™	Short	20 pcs/box
3692908-2	Solo Ultra LiftPants™	XXL, XXXL	20 pcs/box



- 1 SOLO HIGHBACK SLING™**  
is ideal for most lifting situations. It can be used for lifting patients between a bed and wheelchair, or to/from the toilet.
- 2 SOLO REPOSHEET™**  
is used to reposition patients who have moved down in bed, or to turn bed-ridden patients in order to wash them and to help prevent pressure injuries.<sup>2</sup>
- 3 SOLO SUPPORTVEST™**  
is used in conjunction with the Liko Sabina™ sit-to-stand lift, to raise a patient from a sitting to standing position.

- 4 SOLO MULTISTRAPS™**  
are used to facilitate examination, care and treatment. They can also be used to turn a patient in bed and support prone positioning.
- 5 SOLO OCTO LIFTSHEET™**  
is intended for horizontal lifting in conjunction with OctoStretch™. Ideal for use in surgical and intensive care departments.
- 6 SOLO COVER FOR ULTRA LIFTPANTS™**  
is a disposable accessory for use with Liko Ultra LiftPants. The cover is placed over the top of the LiftPants, protecting them from becoming soiled.

<sup>2</sup> Chapter 8: Repositioning and Early Mobilization. In: Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline, The International Guideline 2019. European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance.

For more information, please contact your local distributor or Hillrom sales representative.

[hillrom.com](http://hillrom.com)

Solo disposable products are accessories to be used in situations where individual slings for short term use is needed; for example where there is a risk for cross infection between patients.

Class I  
Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

These products are regulated health care products which, where required by applicable regulations, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documentation accompanying medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.