



Welch Allyn Connex[®] Vital Signs Monitor *with the Connex Scoring App*

Capture a custom patient score at the bedside—right from your Connex Vital Signs Monitor.

The Connex Scoring App provides an optional point-of-care calculator that is customized to your facility's protocols. With the Connex Scoring App, clinicians can document and track signs of patient deterioration and make informed patient care decisions quickly.

WelchAllyn[®]





Customize data fields on the Connex Vital Signs Monitor to match your workflow:

- > Tailor the captured data—specific to a facility or unit
- > Eliminate extra charting steps
- > Up to 20 custom data fields—e.g. integer, decimal, list box or free text
- > Send to your Electronic Medical Record (EMR) to help reduce errors and improve efficiency

Support your facility's early warning score protocol:

- > Automate calculations, such as scores, falls assessments or other data captured at the bedside
- > Create your own custom calculations or match national standards
- > Capture additional data needed, such as intakes, outputs or patient consciousness levels
- > Customize your messages based on calculations
- > Color coding of scores allows for easy recognition of concerns
- > Supports single parameter and aggregate scoring systems



You can configure the device to use whatever protocol you're already using—such as MEWS or other bedside calculations.

For more information about the Connex Vital Signs Monitor, please call 1.800.535.6663 or visit www.welchallyn.com/connex