





# REPLACEMENT BP CUFF REFERENCE GUIDE

|                 | CSM and ProBP 3400  |            |            | CVSM and SPOT LXI  |              |               | Atlas, VSM and Spot Devices (not LXI)   |              |               |
|-----------------|---|------------|------------|--|--------------|---------------|---|--------------|---------------|
| Connector       | No Connector (Raw)  |            |            | 2-tube Locking Type  |              |               | 1-tube Screw-Type   |              |               |
|                 |  |            |            |  |              |               |  |              |               |
| Type            | Durable   | Disposable | Disposable | Durable  | Disposable   | Disposable    | Durable   | Disposable   | Disposable    |
| Material        | Nylon   | Soft       | Vinyl      | Nylon  | Soft         | Vinyl         | Nylon   | Soft         | Vinyl         |
| Unit of Measure | Each  | 20 Pack    | 20 Pack    | Each   | 20 Pack      | 20 Pack       | Each  | 20 Pack      | 20 Pack       |
| Sm. Infant      | Reuse-06  | Soft-06    | Vinyl-06   | N/A  | N/A          | N/A           | Reuse-06-1SC  | Soft-06-1SC  | Vinyl-06-1SC  |
| Infant          | Reuse-07  | Soft-07    | Vinyl-07   | Reuse-07-2MQ   | Soft-07-2MQ  | Vinyl-07-2MQ  | Reuse-07-1SC  | Soft-07-1SC  | Vinyl-07-1SC  |
| Sm. Child       | Reuse-08  | Soft-08    | Vinyl-08   | Reuse-08-2MQ   | Soft-08-2MQ  | Vinyl-08-2MQ  | Reuse-08-1SC  | Soft-08-1SC  | Vinyl-08-1SC  |
| Child           | Reuse-09  | Soft-09    | Vinyl-09   | Reuse-09-2MQ   | Soft-09-2MQ  | Vinyl-09-2MQ  | Reuse-09-1SC  | Soft-09-1SC  | Vinyl-09-1SC  |
| Sm. Adult       | Reuse-10  | Soft-10    | Vinyl-10   | Reuse-10-2MQ   | Soft-10-2MQ  | Vinyl-10-2MQ  | Reuse-10-1SC  | Soft-10-1SC  | Vinyl-10-1SC  |
| Adult           | Reuse-11  | Soft-11    | Vinyl-11   | Reuse-11-2MQ   | Soft-11-2MQ  | Vinyl-11-2MQ  | Reuse-11-1SC  | Soft-11-1SC  | Vinyl-11-1SC  |
| Adult Long      | Reuse-11L   | Soft-11L   | Vinyl-11L  | Reuse-11L-2MQ  | Soft-11L-2MQ | Vinyl-11L-2MQ | Reuse-11L-1SC   | Soft-11L-1SC | Vinyl-11L-1SC |
| Lg. Adult       | Reuse-12  | Soft-12    | Vinyl-12   | Reuse-12-2MQ   | Soft-12-2MQ  | Vinyl-12-2MQ  | Reuse-12-1SC  | Soft-12-1SC  | Vinyl-12-1SC  |
| Lg. Adult Long  | Reuse-12L   | Soft-12L   | Vinyl-12L  | Reuse-12L-2MQ  | Soft-12L-2MQ | Vinyl-12L-2MQ | Reuse-12L-1SC   | Soft-12L-1SC | Vinyl-12L-1SC |
| Thigh           | Reuse-13  | Soft-13    | Vinyl-13   | Reuse-13-2MQ   | Soft-13-2MQ  | Vinyl-13-2MQ  | Reuse-13-1SC  | Soft-13-1SC  | Vinyl-13-1SC  |

|                 | Propaq and PIC Monitors   |              |               | Philips, Hewlett Packard (HP), Spacelabs, Datascope and Mindray Monitors |              |               | Critikon Dinamap Monitors |              |               |
|-----------------|---|--------------|---------------|--|--------------|---------------|---------------------------|--------------|---------------|
| Connector       | 1-tube Locking Type   |              |               | 1-tube Bayonet (aka HP, Rectus)  |              |               | 2-tube Screw-Type         |              |               |
|                 |  |              |               |  |              |               |                           |              |               |
| Type            | Durable   | Disposable   | Disposable    | Durable  | Disposable   | Disposable    | Durable                   | Disposable   | Disposable    |
| Material        | Nylon   | Soft         | Vinyl         | Nylon  | Soft         | Vinyl         | Nylon                     | Soft         | Vinyl         |
| Unit of Measure | Each  | 20 Pack      | 20 Pack       | Each   | 20 Pack      | 20 Pack       | Each                      | 20 Pack      | 20 Pack       |
| Sm. Infant      | Reuse-06-1MQ  | Soft-06-1MQ  | Vinyl-06-1MQ  | Reuse-06-1HP   | Soft-06-1HP  | Vinyl-06-1HP  | N/A                       | N/A          | N/A           |
| Infant          | Reuse-07-1MQ  | Soft-07-1MQ  | Vinyl-07-1MQ  | Reuse-07-1HP   | Soft-07-1HP  | Vinyl-07-1HP  | Reuse-07-2SC              | Soft-07-2SC  | Vinyl-07-2SC  |
| Sm. Child       | Reuse-08-1MQ  | Soft-08-1MQ  | Vinyl-08-1MQ  | Reuse-08-1HP   | Soft-08-1HP  | Vinyl-08-1HP  | Reuse-08-2SC              | Soft-08-2SC  | Vinyl-08-2SC  |
| Child           | Reuse-09-1MQ  | Soft-09-1MQ  | Vinyl-09-1MQ  | Reuse-09-1HP   | Soft-09-1HP  | Vinyl-09-1HP  | Reuse-09-2SC              | Soft-09-2SC  | Vinyl-09-2SC  |
| Sm. Adult       | Reuse-10-1MQ  | Soft-10-1MQ  | Vinyl-10-1MQ  | Reuse-10-1HP   | Soft-10-1HP  | Vinyl-10-1HP  | Reuse-10-2SC              | Soft-10-2SC  | Vinyl-10-2SC  |
| Adult           | Reuse-11-1MQ  | Soft-11-1MQ  | Vinyl-11-1MQ  | Reuse-11-1HP   | Soft-11-1HP  | Vinyl-11-1HP  | Reuse-11-2SC              | Soft-11-2SC  | Vinyl-11-2SC  |
| Adult Long      | Reuse-11L-1MQ   | Soft-11L-1MQ | Vinyl-11L-1MQ | Reuse-11L-1HP  | Soft-11L-1HP | Vinyl-11L-1HP | Reuse-11L-2SC             | Soft-11L-2SC | Vinyl-11L-2SC |
| Lg. Adult       | Reuse-12-1MQ  | Soft-12-1MQ  | Vinyl-12-1MQ  | Reuse-12-1HP   | Soft-12-1HP  | Vinyl-12-1HP  | Reuse-12-2SC              | Soft-12-2SC  | Vinyl-12-2SC  |
| Lg. Adult Long  | Reuse-12L-1MQ   | Soft-12L-1MQ | Vinyl-12L-1MQ | Reuse-12L-1HP  | Soft-12L-1HP | Vinyl-12L-1HP | Reuse-12L-2SC             | Soft-12L-2SC | Vinyl-12L-2SC |
| Thigh           | Reuse-13-1MQ  | Soft-13-1MQ  | Vinyl-13-1MQ  | Reuse-13-1HP   | Soft-13-1HP  | Vinyl-13-1HP  | Reuse-13-2SC              | Soft-13-2SC  | Vinyl-13-2SC  |

|                 | GE Marquette Monitors |              |               |
|-----------------|-----------------------|--------------|---------------|
| Connector       | 2-tube Locking-Type   |              |               |
|                 |                       |              |               |
| Type            | Durable               | Disposable   | Disposable    |
| Material        | Nylon                 | Soft         | Vinyl         |
| Unit of Measure | Each                  | 20 Pack      | 20 Pack       |
| Sm. Infant      | N/A                   | N/A          | N/A           |
| Infant          | Reuse-07-2MQ          | Soft-07-2MQ  | Vinyl-07-2MQ  |
| Sm. Child       | Reuse-08-2MQ          | Soft-08-2MQ  | Vinyl-08-2MQ  |
| Child           | Reuse-09-2MQ          | Soft-09-2MQ  | Vinyl-09-2MQ  |
| Sm. Adult       | Reuse-10-2MQ          | Soft-10-2MQ  | Vinyl-10-2MQ  |
| Adult           | Reuse-11-2MQ          | Soft-11-2MQ  | Vinyl-11-2MQ  |
| Adult Long      | Reuse-11L-2MQ         | Soft-11L-2MQ | Vinyl-11L-2MQ |
| Lg. Adult       | Reuse-12-2MQ          | Soft-12-2MQ  | Vinyl-12-2MQ  |
| Lg. Adult Long  | Reuse-12L-2MQ         | Soft-12L-2MQ | Vinyl-12L-2MQ |
| Thigh           | Reuse-13-2MQ          | Soft-13-2MQ  | Vinyl-13-2MQ  |

