



NOW AVAILABLE
In Four Convenient Sizes

Welch Allyn KleenSpec[®] Single Use LED Vaginal Specula

Integrated light source that is ready to use out of the package,
simply pull switch when ready.



Welch Allyn

KleenSpec Single Use LED Vaginal Specula

NOW AVAILABLE
In Four Convenient Sizes



LED Technology

- > Provides continuous output for 30 minutes
- > Emits crisp white light with no blue or yellow hues for true visualization of tissue color
- > Produces uniform light output, with no hot spots, to help provide the visualization needed for improved patient exams



- > LED light source is sealed within handle to reduce patient risk and placement provides unobstructed view of exam area
- > Smooth molded edges designed to provide maximum patient comfort
- > 100% acrylic, as opposed to an acrylic-styrene mix, is 50% stronger minimizing the chance of cracking or breaking
- > Ready to use out of the package, simply pull switch when ready
- > Look for color-coded Welch Allyn handle colors



- > Lithium-primary battery provides longer shelf life and withstands extreme temperature variations better than alkaline
- > Option to remove lithium-primary battery and LED for disposal or recycling
- > Designed for quality and durability from a trusted manufacturer and made in the U.S.A.



590XS-LED	KleenSpec Single Use LED Vaginal Specula, Extra Small (24/Box, 4 Boxes/Case)
59000-LED	KleenSpec Single Use LED Vaginal Specula, Small (24/Box, 4 Boxes/Case)
59001-LED	KleenSpec Single Use LED Vaginal Specula, Medium (24/Box, 4 Boxes/Case)
59004-LED	KleenSpec Single Use LED Vaginal Specula, Large (14/Box, 4 Boxes/Case)

For more information, please contact your local Welch Allyn representative or call **1.800.535.6663**.

WWW.WELCHALLYN.COM

© 2015 Welch Allyn MC12610

WelchAllyn[®]