

Welch Allyn

PanOptic™ Ophthalmoscope and iExaminer™ System

A mobile digital imaging solution with superior, panoramic views of the fundus



Capture, Store, Send, Retrieve

- iExaminer adapter and PanOptic Ophthalmoscope work together through your iPhone® and the iExaminer App
- Rapid capture of fundus images supports easy, cost-effective eyeground documentation
- Facilitates patient education and compliance
- Compatible with iPhone 6, 6s and 6 Plus models

The Welch Allyn iExaminer App is a required component to capture, store, send and retrieve fundus images. Free and Pro versions are available for Apple iPhone 6 and 6s.




To learn more, visit welchallyn.com/iexaminer.

The Welch Allyn PanOptic Ophthalmoscope is a breakthrough concept in ophthalmoscopy.

Axial PointSource™ optics support fast, easy entry into small, undilated pupils:

- 5X greater field of view versus standard ophthalmoscopes
- **26% increase** in magnification, offering superior views and retinal detail




Patient Eyecup:

- Establishes and maintains proper viewing distance
- Serves as a pivot point for panning around the retina
- Supports orientation and screens out ambient light

Additional features:

- Focus range of -20 to +20 diopters
- Available with Halogen HPX® or SureColor™ LED light source
- Compatible with all Welch Allyn 3.5 V power sources for the desk and wall

Apertures and Filters:



Ordering Information

11810	PanOptic Ophthalmoscope head only (includes: micro, small, large spot; slit and half-moon apertures; red-free filter) and Halogen HPX lamp
11820	PanOptic Ophthalmoscope head only (includes: micro, small, large spot; slit aperture; red-free and cobalt blue filter; magnification lens) and Halogen HPX lamp
11820-L	PanOptic Ophthalmoscope head only (includes: micro, small, large spot; slit aperture; red-free and cobalt blue filter; magnification lens) and SureColor LED lamp
11840	iExaminer adapter for iPhone 4 and 4s
11840-A6	iExaminer adapter for iPhone 6 and 6s
11840-A6P	iExaminer adapter for iPhone 6 Plus and 6s Plus
11842	PanOptic iExaminer kit for 4 and 4s
11842-A6	PanOptic iExaminer kit for 6 and 6s
11842-A6P	PanOptic iExaminer kit for 6 Plus and 6s Plus

Accessories

11870	Reusable patient eyecups (package of 5)
11875	Corneal viewing lens (for use with 11820 and 11820-L)
03800-U	03800-U 3.5 V Halogen HPX replacement lamp
03800-LED	SureColor LED replacement lamp

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM