



## Welch Allyn CP 50™ Electrocardiograph



Simple. Affordable. Compact.

The practical CP 50 ECG—designed to improve office efficiencies.


- > Offered at an affordable price to help maximize your budget
- > Easy-to-use, with one-button operation, makes performing ECG tests efficient
- > High-resolution thermal printer provides quick outputs
- > Full-color, touch-screen display to help improve staff satisfaction
- > Battery operation lets you easily transport your device to the patient
- > Optional Mobile Stand and Carrying Case are also available
- > Optional adult and pediatric MEANS interpretation offers a second opinion



EASY TO USE

### Touch-Screen Display


Facilitates one-button operation, guiding you through each test quickly and easily.



EASY TO FOLLOW

### On-Screen Instructions

Saves valuable time and helps ensure accurate data capture by highlighting poor signals and connections.



EASY TO READ

### Programmable Report Formats

Prints reports in the format that works best for you with a high-resolution thermal printer.

**WelchAllyn®**


Advancing Frontline Care™

# Improve your practice workflow and efficiency with the Welch Allyn CP 50™ ECG



## Lightweight and portable

- One-button operation facilitates quick setup and operation
- Optional carrying case
- Optional mobile stand
- Weighs just over 4.4 lbs



Welch Allyn CardioPerfect™ Workstation

## Bring your cardiopulmonary data together with the Welch Allyn CardioPerfect® Workstation

- Manage patient information and tests performed with the full suite of Welch Allyn cardiopulmonary devices
- Get paperless efficiency
- Store, access, review, compare and edit diagnostic test data from any computer on your network
- Retrieve records with the touch of a key, and review and edit right on the screen
- Access information from your exam rooms, your office or at home



CP 50 Rolling Cart, Part # 406815

## Ordering/Parts information for the CP 50

### Device configurations:

- CP50-1EN1 CP 50—Non-interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.
- CP50A-1EN1 CP 50—Interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.
- CP50P-1EN1 CP 50 Plus—Non-interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.
- CP50AP-1EN1 CP 50 Plus—Interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.

### Accessories:

- |            |   |        |  |
|------------|---|--------|--|
| 715006     | ECG Multi-function Electrode Adapter      | 406020 | CP 50 Printer Paper, Roll 4 rolls/case |
| 45008-0000 | Resting Tab Electrodes (box of 1000)      | 715107 | CP 50 Carrying Case                    |
| 719653     | 10-lead Patient Cable, AHA, Banana, CP 50 | 406815 | CP 50 Cart                             |
| 406021     | CP 50 Printer Paper Z fold 4 packs/case   |        |  |

### Partners in Care<sup>SM</sup> Service Options:

- S1-CP50 CP 50, Comprehensive Partner Program
- S1-CP50-2 CP 50, Comprehensive Partner Program 2YR

### Warranty information

- CP 50 device—1 year
- Patient Cable /Battery—90 days

## Flexible configuration options to meet the needs of your practice

Choose the CP 50 that's best for you



### Welch Allyn CP 50™

- Reliable touch-screen display
- Optional MEANS interpretation
- High-resolution thermal printer
- Internal storage

### Welch Allyn CP 50™ Plus with Connectivity

Includes all of the features of the CP 50, plus:

- Flexible connectivity options
- External printing capability
- Programmable patient entry fields
- Administrative controls
- Store ECGs to a USB memory stick as PDFs

**WelchAllyn®**

Advancing Frontline Care™

Contact your local Welch Allyn or distributor representative for a full list of configurations, or visit [www.welchallyn.com/CP50](http://www.welchallyn.com/CP50) today for more information.

Welch Allyn Corporate Headquarters

4341 State Street Road, Skaneateles Falls, NY 13153-0220 USA

Telephone: +1.800.535.6663

[www.welchallyn.com/CP50](http://www.welchallyn.com/CP50)

© 2013 Welch Allyn SM4098 Rev A