

Welch Allyn TM55 and TM65 Treadmills

Specifications

Primary Benefits:

- **Ease of use** — The treadmill is easy to operate and whisper quiet at all speeds.
- **Control** — The Ultradrive system is designed for smooth power and accuracy even at the slowest speeds and the heaviest weights.
- **Integration** — TM55 and TM65 treadmills are compatible with Q-Stress® and XScribe™ stress testing systems.



The TM55 and TM65 treadmills provide a comfortable and quiet treadmill experience.

Clinicians around the world rely on our treadmills. Your patients will appreciate the solid and secure feeling the TM55 and TM65 treadmills provide, while you will appreciate how easy they are to use.

- Our easy-to-read cardiac stress systems make it simple for you to start, control and stop the treadmill.
- Clear displays let you quickly find test duration, speed and grade.
- The stress system allows you to bring your patients to a comfortable stop, gradually reducing speed and lowering the grade to zero.
- Optional rapid deceleration capability is available for the times you need to slow and stop the treadmill quickly, such as stress echo peak exercise imaging.

The Ultradrive system is the power behind the TM55 and TM65 treadmills.

Ultradrive is a digitally controlled AC-drive system that provides exceptional performance, accuracy and reliability.

- Treadmills deliver high torque for accurate speed control, even at lower speeds.
- Patients enjoy a natural walking/running sensation with no lag when the foot strikes the belt.

Welch Allyn TM55 and TM65 Treadmills

Technical Specifications

Dimensions	<p>Overall: (H x L x W)</p> <p>TM55 — 46.5 in x 80.2 in x 29.9 in (118 cm x 204 cm x 75.9 cm)</p> <p>TM65 — 46.5 in x 90.2 in x 29.9 in (118 cm x 229 cm x 75.9 cm)</p> <p>Running Surface: (H x L x W)</p> <p>TM55 — 6.3 in x 55 in x 20 in (16 cm x 140 cm x 51 cm)</p> <p>TM65 — 6.3 in x 65 in x 20 in (16 cm x 165 cm x 51 cm)</p> <p>Optional Side Handrails: (H x L)</p> <p>Short — 24 in x 40 in (61 cm x 102 cm)</p> <p>Long — 44 in x 40 in (112 cm x 102 cm)</p>
Weight	<p>Treadmill:</p> <p>TM55 — 352 lb (160 kg)</p> <p>TM65 — 375 lb (170 kg)</p> <p>Shipping:</p> <p>TM55 — 500 lb (227 kg)</p> <p>TM65 — 525 lb (238 kg)</p>
User Weight Capacity	500 lb (227 kg)
Speed Range	0.8 mph to 9.6 mph (1.3 km/h to 15.4 km/h) in 0.1 mph increments
Incline Range	0% to 25%, in 0.5% increments
Warranty	13 months parts and labor
Power	<p>Low Voltage — 100/120 VAC @ 20 A; 50/60 Hz; NEMA 5-20R non-locking receptacle; Single-phase dedicated circuit</p> <p>High Voltage — 200/230 VAC @ 10 A; 50/60 Hz; NEMA 6-15R non-locking receptacle; Single-phase dedicated circuit</p>

Specifications subject to change without notice.

For more information, contact your local Welch Allyn representative or visit www.welchallyn.com.



Clinical excellence. Connected solutions. Continuous innovation.
Welch Allyn Cardiology is proud to be powered by Mortara.

